

pGL3 Luciferase Reporter Vectors

Technical Manual No. 033

INSTRUCTIONS FOR USE OF PRODUCTS E1741, E1751, E1761, E1771.

PLEASE DISCARD PREVIOUS VERSIONS.

All Technical Literature is Available on the Internet at www.promega.com

Please visit the web to verify that you are using the most current version of this Technical Manual.

I. Description	1
II. pGL3 Vector Maps and Sequence Reference Points	2
III. Product Components	6
IV. Cloning Methods.....	6
A. Cloning Strategies	6
B. Preparation of pGL3 Vectors and Insert DNA for Cloning	6
C. Transformation Protocols for pGL3 Vectors	7
D. Isolation of Plasmid DNA.....	7
V. Transfection of Mammalian Cells	8
VI. Assay of Luciferase Activity	9
VII. Generation of Nested Deletions	10
VIII. Generation of Single-Stranded DNA and Site-Specific Mutations	10
A. Production of Single-Stranded DNA	10
B. Generation of Site-Specific Mutations	11
IX. Sequencing of Luciferase Reporter Vectors	11
X. Appendix	12
A. Common Structural Elements of the pGL3 Luciferase Reporter Vectors	12
B. Advantages of the pGL3 Vectors	13
C. Description of Reporter Vector Changes	13
D. Distinguishing Features of the pGL3 Luciferase Reporter Vectors.....	15
E. Mapping Genetic Elements Located Within DNA Fragments	15
F. Composition of Buffers and Solutions	16
G. References	17
H. pGL3-Basic Vector Restriction Sites and Sequence.....	19
I. pGL3-Enhancer Vector Restriction Sites and Sequence	24
J. pGL3-Promoter Vector Restriction Sites and Sequence.....	29
K. pGL3-Control Vector Restriction Sites and Sequence.....	34
<i>Experienced User's Protocol</i>	40

I. Description

The pGL3 Luciferase Reporter Vectors^(a,b) provide a basis for the quantitative analysis of factors that potentially regulate mammalian gene expression. These factors may be cis-acting, such as promoters and enhancers, or trans-acting, such as various DNA-binding factors. The backbone of the pGL2 Luciferase Reporter Vectors^(b) was redesigned for the pGL3 Vectors for increased expression, and

contains a modified coding region for firefly (*Photinus pyralis*) luciferase that has been optimized for monitoring transcriptional activity in transfected eukaryotic cells. The assay of this genetic reporter is rapid, sensitive and quantitative. In addition, the Luciferase Reporter Vectors contain numerous features aiding in the structural characterization of the putative regulatory sequences under investigation.

II. pGL3 Vector Maps and Sequence Reference Points




Figure 1. pGL3-Basic Vector circle map. Additional description: *luc+*, cDNA encoding the modified firefly luciferase; *Amp^r*, gene conferring ampicillin resistance in *E. coli*; f1 ori, origin of replication derived from filamentous phage; ori, origin of replication in *E. coli*. Arrows within *luc+* and the *Amp^r* gene indicate the direction of transcription; the arrow in the f1 ori indicates the direction of ssDNA strand synthesis.

pGL3-Basic Vector Sequence Reference Points:

SV40 Promoter	(none)
SV40 Enhancer	(none)
Multiple cloning region	1–58
Luciferase gene (<i>luc+</i>)	88–1740
GLprimer2 binding site	89–111
SV40 late poly(A) signal	1772–1993
RVprimer4 binding site	2080–2061
Co/E 1-derived plasmid replication origin	2318
β-lactamase gene (<i>Amp^r</i>)	3080–3940
f1 origin	4072–4527
Synthetic poly(A) signal	4658–4811
RVprimer3 binding site	4760–4779




Figure 2. The pGL3-Enhancer Vector circle map. Additional description: *luc+*, cDNA encoding the modified firefly luciferase; *Amp^r*, gene conferring ampicillin resistance in *E. coli*; f1 ori, origin of replication derived from filamentous phage; ori, origin of plasmid replication in *E. coli*. Arrows within *luc+* and the *Amp^r* gene indicate the direction of transcription; the arrow in f1 ori indicates the direction of ssDNA strand synthesis.

pGL3-Enhancer Vector Sequence Reference Points:

SV40 Promoter	(none)
Multiple cloning region	1–58
Luciferase gene (<i>luc+</i>)	88–1740
GLprimer2 binding site	89–111
SV40 late poly(A) signal	1772–1993
SV40 Enhancer	2013–2249
RVprimer4 binding site	2307–2326
<i>ColE1</i> -derived plasmid replication origin	2564
β -lactamase gene (<i>Amp^r</i>)	3326–4186
f1 origin	4318–4773
Synthetic poly(A) signal	4904–5057
RVprimer3 binding site	5006–5025




Figure 3. The pGL3-Promoter Vector circle map. Additional description: *luc+*, cDNA encoding the modified firefly luciferase; *Amp^r*, gene conferring ampicillin resistance in *E. coli*; f1 ori, origin of replication derived from filamentous phage; ori, origin of replication in *E. coli*. Arrows within *luc+* and the *Amp^r* gene indicate the direction of transcription; the arrow in f1 ori indicates the direction of ssDNA strand synthesis.

pGL3-Promoter Vector Sequence Reference Points:

Multiple cloning region	1–41
SV40 Promoter	48–250
Luciferase gene (<i>luc+</i>)	280–1932
GLprimer2 binding site	281–303
SV40 Enhancer	(none)
SV40 late poly(A) signal	1964–2185
RVprimer4 binding site	2253–2272
<i>Co/E</i> 1-derived plasmid replication origin	2510
β -lactamase gene (<i>Amp^r</i>)	3272–4132
f1 origin	4264–4719
Synthetic poly(A) signal	4850–5003
RVprimer3 binding site	4952–4971




Figure 4. pGL3-Control Vector circle map. Additional description: *luc+*, cDNA encoding the modified firefly luciferase; *Amp^r*, gene conferring ampicillin resistance in *E. coli*; f1 ori, origin of replication derived from filamentous phage; ori, origin of replication in *E. coli*. Arrows within *luc+* and the *Amp^r* gene indicate the direction of transcription; the arrow in f1 ori indicates the direction of ssDNA strand synthesis.

pGL3-Control Vector Sequence Reference Points:

Multiple cloning region	1-41
SV40 Promoter	48-250
Luciferase gene (<i>luc+</i>)	280-1932
GLprimer2 binding site	281-303
SV40 late poly(A) signal	1964-2185
SV40 Enhancer	2205-2441
RVprimer4 binding site	2499-2518
<i>ColE1</i> -derived plasmid replication origin	2756
β -lactamase gene (<i>Amp^r</i>)	3518-4378
f1 origin	4510-4965
Synthetic poly(A) signal	5096-5249
RVprimer3 binding site	5198-5217




Figure 5. pGL3 Vector multi-cloning regions. The upstream and downstream cloning sites and the location of the sequencing primers, GLprimer2, RVprimer3 and RVprimer4 are shown. The large primer arrows indicate the direction of sequencing. The positions of the promoter (in the pGL3-Promoter and pGL3-Control Vectors) and the enhancer (in the pGL3-Enhancer and pGL3-Control Vectors) are shown as insertions into the sequence of the pGL3-Basic Vector. (Note that the promoter replaces four bases [AAGT] of the pGL3-Basic Vector.) The sequence shown is of the DNA strand generated from the f1 ori.

III. Product Components

Product	Size	Cat.#
pGL3-Control Vector(a,b)	20µg	E1741
pGL3-Basic Vector(a,b)	20µg	E1751
pGL3-Promoter Vector(a,b)	20µg	E1761
pGL3-Enhancer Vector(a,b)	20µg	E1771

Vectors are supplied with a glycerol stock of bacterial strain JM109 cells. The JM109 cells do not contain the vector and are not competent cells. Information on related products, including the Luciferase Assay System, is provided in Sections III–VIII.

Product	Size	Cat.#
GLprimer2 (counter clockwise)	2µg	E1661
RVprimer3 (clockwise)	2µg	E4481
RVprimer4 (counter clockwise)	2µg	E4491

Storage Conditions: Store the pGL3 Luciferase Reporter Vectors at –20°C and the glycerol stock of JM109 cells at –70°C.

IV. Cloning Methods

A. Cloning Strategies

The restriction sites for *Xho* I and *Sal* I have compatible ends, as do *Bgl* II and *Bam*H I. Therefore, cloning into the *Xho* I or *Bgl* II sites upstream of *luc+*, or the downstream *Sal* I or *Bam*H I sites, allows for easy interchange of DNA inserts between upstream and downstream positions relative to the luciferase reporter gene. Thus, positional effects of a putative genetic element may be readily tested. Cloning fragments into a single site will generally yield both possible orientations relative to the reporter gene, making these effects also readily examinable.

The other upstream restriction sites may be used for cloning. However, note that some of the sites are required for generation of nested deletions (see Section VII). Specifically, the *Kpn* I or *Sac* I site is needed to generate a 3'-overhang upstream of the insert.

B. Preparation of pGL3 Vectors and Insert DNA for Cloning

The fragment and vector DNA should be digested with restriction enzymes that will generate compatible ends for cloning. In some cases, the ends of the DNA fragment may require modification, either by using synthetic linkers, by a PCR^(c) amplification using primers containing sites for appropriate restriction enzymes, or by filling in the restriction site overhangs. It may be advantageous to treat the vector DNA with Calf Intestinal Alkaline Phosphatase (CIAP; Cat.# M2825) or Shrimp Alkaline Phosphatase (Cat.# M8201) to remove 5' phosphate groups, thus preventing reclosure of the vector on itself without an insert. Sufficient DNA should be prepared to perform control reactions for digestion, ligation and transformation steps.

To ensure capture of the correct insert DNA, the desired restriction fragment can be purified by electrophoresis on an acrylamide or agarose gel and then recovered from the gel by one of several methods, such as using the *Wizard*[®] PCR Preps DNA Purification System^(d) Technical Bulletin, #TB118. Alternatively, nonfractionated restriction fragments can be cloned into the target plasmid, and the desired recombinant can then be identified by gel electrophoresis of plasmid DNA.

Protocols for restriction digestion, alkaline phosphatase treatment, linker ligation and transformation of competent cells can be found in Promega's *Protocols and Applications Guide* (1) or in *Molecular Cloning, A Laboratory Manual* (2).

C. Transformation Protocols for pGL3 Vectors

Because the Luciferase Reporter Vectors are supplied as modified DNA, *E. coli* hosts may be either restriction + or restriction -. The Luciferase Reporter Vectors are supplied with JM109 bacterial cells (*endA1*, *recA1*, *gyrA96*, *thi*, *hsd R17*, (*rK-*, *mK+*), *relA1*, *supE44*, Δ (*lac-proAB*), [*F'*, *traD36*, *proAB*, *lac ρ Z Δ M15*]). The use of a *recA* host such as JM109 is preferred because this prevents undesirable recombination between the insert and the host chromosomal DNA. A strain that has an *F'* episome is required for ssDNA production.

Grow JM109 on minimal plates (M-9) supplemented with 1.0mM thiamine-HCl prior to preparation of competent cells and transformation. This selects for the presence of the *F'* episome.

D. Isolation of Plasmid DNA

The *Wizard*[®] Plus *SV Minipreps DNA Purification System*^(e) *Technical Bulletin*, #TB225, or the *Wizard*[®] Plus *Midipreps DNA Purification System*^(f) *Technical Bulletin*, #TB173, may be used for small-scale preparation of plasmid DNA for screening clones. Large-scale DNA preparations can be made for sequencing or restriction digestion using the *Wizard*[®] Plus *Maxipreps DNA Purification System*^(f) *Technical Bulletin*, #TB139, or the *Wizard*[®] Plus *Megapreps DNA Purification System*^(f) *Technical Bulletin*, #TB140. DNA suitable for transfection may be purified using a modification of the *Wizard*[®] Maxipreps protocol (3) or by using the *Wizard*[®] *PureFection Plasmid DNA Purification System*^(g) *Technical Bulletin*, #TB259, or by CsCl gradient preparation (4).

The following protocol allows for the rapid isolation of large quantities of plasmid DNA without the need for column purification or banding in CsCl gradients (4). The procedure takes advantage of the rapid alkaline denaturation of plasmid and chromosomal DNA and the selective renaturation of plasmid DNA following neutralization of the solution. The polyethylene glycol (PEG) precipitation step is included to help remove contaminants that could interfere with restriction digestions, sequencing procedures or transfection of mammalian cells.

The volume of the culture used may be adjusted depending upon the amount of DNA required for subsequent manipulations. The Luciferase Reporter Vectors are high copy number plasmids (200–300 copies per cell). Using the protocol described below, 150–500 μ g of DNA can be obtained from a 250ml culture.

Materials to Be Supplied by the User

(Solution compositions are provided in Section X.F.)

- LB medium
- lysis buffer for plasmid preps
- TE-saturated phenol:chloroform: isoamyl alcohol
- TE buffer
- 13% (w/v) polyethylene glycol (M.W. 6,000–8,000) in water
- potassium acetate solution (pH 4.8)

Note: An Experienced User's Protocol can be found at the end of this Technical Manual.

1. Prepare 250ml of culture by incubating overnight in LB medium containing 100µg/ml ampicillin.
2. Centrifuge the cells at 5,000 × *g* for 15 minutes at 4°C. Remove and discard the supernatant.
3. Resuspend the cells in 6ml of freshly prepared ice-cold lysis buffer by careful pipetting with a 10ml pipette. Incubate in ice water for 20 minutes.
4. Add 12ml of 0.1N NaOH, 1% SDS (prepared fresh). Mix carefully and thoroughly by inversion. Do not vortex.
5. Add 7.5ml of potassium acetate solution (pH 4.8). Mix carefully by inversion and incubate in ice water for 10 minutes.
6. Centrifuge at 12,000 × *g* for 15 minutes. Transfer the supernatant to a fresh tube, avoiding the white precipitate. Add 50µl of RNase A (1mg/ml stock) to the supernatant. Incubate for 20 minutes at 37°C.
7. Extract with 1 volume of TE-saturated phenol:chloroform:isoamyl alcohol. Centrifuge at 12,000 × *g* for 10 minutes.
8. Save the upper, aqueous phase and repeat the TE-saturated phenol:chloroform:isoamyl alcohol extraction as described in Step 7 above.
9. Extract with one volume of chloroform:isoamyl alcohol (24:1) by vortexing for 1 minute. Centrifuge at 12,000 × *g* for 10 minutes.
10. Transfer the upper, aqueous phase to a fresh tube and add 2 volumes of 100% ethanol. Centrifuge at 12,000 × *g* for 20 minutes.
11. **Optional:** Dissolve the pellet in 1.6ml of water. Add 0.4ml of 4M NaCl and mix. Add 2ml of 13% (w/v) polyethylene glycol (PEG, M.W. 6,000–8,000) and mix. Incubate in ice water for 60 minutes. PEG is used to separate small nucleotides from plasmid DNA. Centrifuge at 12,000 × *g* for 10 minutes.
12. Remove the supernatant and wash the pellet with 70% ethanol. Centrifuge at 12,000 × *g* for 5 minutes.
13. Dry the pellet under vacuum. Dissolve the pellet in water or TE buffer (100–500µl).

V. Transfection of Mammalian Cells

Transfection of DNA into eukaryotic cells may be mediated by cationic lipid compounds (5), calcium phosphate (6,7), DEAE-dextran (6,8), or electroporation (7). Transfection systems based on cationic lipids (e.g., Transfectam[®] Reagent^(h) [Cat.# E1232], TransFast[™] Reagent⁽ⁱ⁾ [Cat.# E2431], or Tfx[™] Reagents^(j) [Cat.# E1811, E2381, E2391]), calcium phosphate and DEAE-dextran are available from Promega. For information on the Transfectam[®] protocol, please request the *Transfectam[®] Reagent Technical Bulletin*, #TB116. For information regarding use of the TransFast[™] Transfection Reagent, request Technical Bulletin #TB260. Protocols for the use of Tfx[™] Reagents are included in Technical Bulletin #TB216. For transfection procedures using calcium phosphate or DEAE-dextran, please request the *ProFection[®] Mammalian Transfection System Technical Manual*, #TM012 (Cat.# E1200, E1210).

VI. Assay of Luciferase Activity

Experimental strategies using firefly luciferase may involve the analysis of a few samples per day or as many as several thousand samples per hour, and equipment used to measure luminescence may vary from inexpensive, single-sample luminometers to high-end CCD luminometers. To support this wide range of applications, Promega has developed three luciferase assays with different, but complementary, characteristics: Luciferase Assay System^(k) (Cat.# E1500), Bright-Glo[™] Luciferase Assay System^(k) (Cat.# E2610), and Steady-Glo[®] Luciferase Assay System^(k) (Cat.# E2510). Reagent choice depends on weighing the relative importance of experimental format, assay sensitivity, and luminescence duration.

Table 1. Characteristics of Promega's Luciferase Assay Reagents.

	Bright-Glo[™] Reagent	Steady-Glo[®] Reagent	Luciferase Assay Reagent
Format	NH or H	NH or H	NH
Process	continuous	batch	bench scale
Number of Steps	1	1	4
Sensitivity	highest	lower	higher
Signal Half-Life	~30 minutes	~5 hours	~12 minutes
Precision	High	High	High
Cell Lysis Time	~2 minutes maximum	~5 minutes maximum	NA
Reagent Prep Time	<30 seconds	<30 seconds	Up to 40 minutes

NH = nonhomogeneous; H = homogeneous; NA = not applicable

The Luciferase Assay System has long been the standard reagent for routine laboratory analysis. Before using this reagent, cells from which the luciferase is to be measured must be washed and lysed. This reagent was optimized for high sensitivity in nonhomogeneous, single-sample measurements. The Luciferase Assay System requires a luminometer fitted with injectors to efficiently measure luminescence in 96-well plates.

The Bright-Glo[™] and Steady-Glo[®] Reagents were developed to perform assay reactions within multiwell plates and in the presence of complete cell culture medium: no cell preparation steps such as washing or lysing are required before the luminescence reaction is initiated. Both of these are single-step reagents, requiring only addition of the reagent before measuring luminescence. This makes them ideal reagents for efficient and precise quantitation in 96-, 384- and 1536-well plates.

The Bright-Glo[™] and Steady-Glo[®] Reagents are complementary in their characteristics based on the inverse relationship between luminescence duration and assay sensitivity (9). Generally as the half-life of the luminescence increases, assay sensitivity decreases. The Steady-Glo[®] Reagent provides very long luminescence duration (changing only about 10% per hour); however, to achieve this long luminescence duration, the assay sensitivity must be reduced. This reagent was designed for experimental designs in which many microplates are processed as a batch.

In contrast, the Bright-Glo™ Reagent provides high assay sensitivity with lower luminescence duration (<10% decrease per 5 minutes). This reagent is designed for general research applications and for experimental designs using robotics for continuous sample processing. Furthermore, as a result of increased sample capacity, the Bright-Glo™ Reagent provides greater assay sensitivity than the Luciferase Assay Reagent in most applications (9).

The Luciferase Assay System, Bright-Glo™ Reagent and Steady-Glo® Reagent provide the highest standards in assay quantitation, sensitivity and convenience. Since these reagents are based on the same underlying design principles, different reagents can be used as experimental needs change. For more information request the *Luciferase Assay System Technical Bulletin*, #TB281, the *Steady-Glo® Luciferase Assay System Technical Manual*, #TM051, or the *Bright-Glo™ Luciferase Assay System Technical Manual*, #TM052.

VII. Generation of Nested Deletions

Unidirectional deletions of any inserted DNA can be made using a procedure developed by Henikoff (10) in which Exonuclease III (Exo III) is used to specifically digest insert DNA from a 5' protruding or blunt-ended restriction site. In the pGL3 Luciferase Reporter Vectors, these 5' overhangs are supplied by digesting the plasmid with *Bgl* II, *Mlu* I, *Nhe* I, *Xho* I or *Xma* I. When the plasmids are cut with *Kpn* I or *Sac* I, which yield 3' overhangs, the Exo III will be unable to digest in the other direction.

The uniform rate of enzyme digestion allows deletions of various lengths to be made simply by removing timed aliquots from the reaction. Given that small deletions (less than 500 bases) are probably desired, we recommend performing the reactions at a lower temperature (between 4–16°C). Samples from the Exo III reaction are removed at timed intervals to tubes containing S1 nuclease, which removes the remaining single-stranded tails. The low pH and the presence of zinc cations in the S1 buffer effectively inhibit further digestion by Exo III. After neutralization and heat inactivation of the S1 nuclease, Klenow DNA polymerase is added to flush the ends, which are then ligated to circularize the deletion-containing vectors. The ligation mixtures are used directly to transform competent cells. Each successive time point yields a collection of subclones containing clustered deletions extending further into the original insert.

For a more detailed protocol, please request the *Erase-a-Base® System Technical Manual*, #TM006.

VIII. Generation of Single-Stranded DNA and Site-Specific Mutations

A. Production of Single-Stranded DNA

To generate single-stranded DNA (ssDNA) from the pGL3 Vectors, bacterial cells containing pGL3 Vectors are infected with an appropriate helper phage. The plasmid then enters the f1 replication mode, and the resulting ssDNA is exported from the cell as an encapsulated phage-like particle. The single-stranded plasmid DNA is purified from the supernatant by simple precipitation and extraction procedures. Promega's *Protocols and Applications Guide* (1) contains protocols for the preparation and analysis of ssDNA suitable for mutagenesis and sequencing (1,11–13).

B. Generation of Site-Specific Mutations

Site-specific mutagenesis, as developed by Hutchinson *et al.* (14), is accomplished by hybridizing to ssDNA a synthetic oligonucleotide that is complementary to the single-stranded template except for a region of mismatch near the center. It is this region that contains the desired nucleotide change or changes. Following hybridization with the single-stranded target DNA, the oligonucleotide is extended with DNA polymerase to create a double-stranded structure. The nick is then sealed with DNA ligase, and the duplex structure is transformed into an *E. coli* host. Theoretically, the yield of mutants using the Hutchinson procedure should be 50% (due to semi-conservative replication). In practice, however, the mutant yield may be much lower, often only a few percent or less. This is assumed to be due to factors such as incomplete *in vitro* polymerization, primer displacement by the DNA polymerase used in the fill-in reaction, and *in vivo* host-directed mismatch repair mechanisms, which favor repair of the nonmethylated, newly synthesized DNA strand. Because of the low mutant yield, methods have been developed to increase the mutation frequency.

Promega's Altered Sites® II *in vitro* Mutagenesis Systems^(l) (Cat.# Q6080, Q6090, Q6210) use antibiotic selection to obtain consistently high mutagenesis frequencies (often >90%) using ssDNA or dsDNA templates. These systems provide a simple, one-day procedure for generation and selection of oligonucleotide-directed mutants and include the ability to perform sequential rounds of mutagenesis without subcloning and to express the mutated gene products *in vivo* or *in vitro*. For further information, please request the *Altered Sites® II in vitro Mutagenesis System Technical Manual*, #TM001.

IX. Sequencing of Luciferase Reporter Vectors

It may be desirable to sequence the DNA inserted into the Luciferase Reporter Vectors. Two examples of such applications are to determine the exact position of generated deletions (see Section VII) and to confirm production of a site-specific mutation (see Section VIII.B). Three primers are available for sequencing the pGL3 Vectors: RVprimer3 (Reporter Vector Primer 3) for sequencing clockwise across the upstream cloning sites, RVprimer4 for sequencing counterclockwise across the *Bam*H I and *Sal*I cloning sites downstream of *luc+*, and GLprimer2 for sequencing counterclockwise upstream of *luc+*.

RVprimer3	5'-CTAGCAAATAGGCTGTCCC-3'
RVprimer4	5'-GACGATAGTCATGCCCCGCG-3'
GLprimer2	5'-CTTTATGTTTTTGGCGTCTTCCA-3'

RVprimer3 is especially useful for identifying positions of nested deletions. Note that all three primers can be used for dsDNA sequencing, but only RVprimer4 and GLprimer2 can also be used for ssDNA sequencing.

There are many methods for DNA sequencing; the most appropriate method will depend on the specific application and on your experience. Thermal cycle sequencing takes advantage of the intrinsic properties of the DNA polymerase isolated from *Thermus aquaticus* (*Taq* DNA polymerase). *Thermus aquaticus* is an extremely thermophilic microorganism whose DNA polymerase shows thermal stability to 95°C (15,16). Promega's *fmoI*^(m,o) and SILVER SEQUENCE™ DNA Sequencing Systems^(m,n,o) include Promega's Sequencing Grade *Taq* DNA Polymerase^(o)

(Cat.# M2031, M2035), which is a modified form that gives superior results on dsDNA templates because it lacks 5'→3' exonuclease activity. Sequencing Grade *Taq* DNA Polymerase produces a uniform band intensity, low background and a high degree of accuracy.

These sequencing systems use a thermocycling apparatus that yields a linear amplification of template DNA, decreasing the amount of template necessary to obtain sequence data. The high temperatures used in this procedure eliminate the need for alkaline denaturation and ethanol precipitation of dsDNA templates. They also increase the stringency of primer hybridization, providing more accurate sequence data, and decrease DNA secondary structure, permitting polymerization through highly structured regions (16). The *fmol*[®] System (Cat.# Q4100) provides the option of using either ³²P or ³⁵S radioactively end-labeled primers or direct incorporation of radioactive label for sequencing reactions. The SILVER SEQUENCE[™] System (Cat.# Q4130), by contrast, uses silver staining detection to eliminate the need for radioactivity in sequencing reactions.

For further information, please request the *fmol*[®] or SILVER SEQUENCE[™] Sequencing System Technical Manuals, #TM024 or #TM023, respectively.

X. Appendix

A. Common Structural Elements of the pGL3 Luciferase Reporter Vectors

Except for the promoters and enhancers, the four pGL3 Luciferase Reporter Vectors are structurally identical. Each plasmid's distinguishing features are summarized in Section X.D. The pGL3 Vectors each contain a high copy number prokaryotic origin of replication for maintenance in *E. coli*, an ampicillin-resistance gene for selection, and a filamentous phage origin of replication (f1 ori) for single-stranded DNA (ssDNA) production. Restriction sites for insertion of DNA fragments are located upstream and downstream of the luciferase gene. Two of the upstream sites (*Xho* I and *Bgl* II) yield cohesive ends compatible with the downstream sites (*Sal* I and *Bam*H I, respectively), allowing the interchange of the DNA insert for rapid analysis of positional effects.




Figure 6. Comparison of luciferase activities expressed in HeLa cells transfected with the pGL2-Control and pGL3-Control Reporter Vectors. The expression level of *Luc+* is dramatically higher with the pGL3-Control Vectors. In repeated experiments with several cell lines, we observed 20- to 100-fold higher luciferase activity from cells transfected with pGL3-Control. Luciferase activity was measured with a Turner Designs luminometer. (Absolute light values and relative expression profiles may vary between different cell types.)




Figure 7. A representative experiment comparing luciferase activities expressed in HeLa cells transfected with the pGL2 and pGL3 Vector series. The increase in luciferase expression observed with these new vectors provides greater sensitivity, while maintaining relatively low background luciferase expression.

B. Advantages of the pGL3 Vectors

The pGL3 Luciferase Reporter Vectors provide significant advances over the pGL2 Reporter Vectors. The pGL3 Reporter Vectors contain a modified firefly luciferase cDNA designated *luc+* and a redesigned vector backbone. These changes increase luciferase expression, improve in vivo vector stability, and provide greater flexibility in performing genetic manipulations. The modified reporter vectors result in luciferase expression levels dramatically higher than those obtained with pGL2 Reporter Vectors (Figure 6), while maintaining relatively low background luciferase expression (Figure 7).

The substantial increase in the expression of luciferase observed with these new vectors provides greater sensitivity. It may now be possible to obtain measurable luciferase expression in cell types that are difficult to transfect or when studying weak promoter elements. Users of the pGL2 and pGL3 Vectors should be aware, however, that absolute light unit values and relative expression profiles vary between different cell types (17). Therefore, it is important to include the appropriate control vectors in all experiments.

C. Description of the Reporter Vector Changes

Modifications were made to both the luciferase gene (*luc+*) and the vector backbone. The modifications that distinguish the *luc+* gene from the native luciferase gene generally fall into four categories: i) the C-terminal tripeptide has been removed to eliminate peroxisome targeting of the expressed protein; ii) codon usage was improved for expression in plant and animal cells; iii) two potential sites of N-glycosylation were removed; and iv) several DNA sequence changes were made to disrupt extended palindromes, remove internal restriction sites, and eliminate consensus sequences recognized by genetic regulatory binding proteins, thus helping to ensure that the reporter gene itself is unaffected by spurious host transcriptional signals. (For a detailed description of the modifications to the *luc+* gene, refer to #TB208 and reference 19.)

Four major modifications were made to the pGL2 vector backbone: i) the SV40 early poly(A) signal has been replaced with the SV40 late poly(A) signal to increase the efficiency of transcription termination and polyadenylation of the luciferase transcripts (19); ii) a synthetic poly(A) and transcriptional pause site (20,21) have been placed upstream of the multiple cloning site to terminate spurious transcription, which may initiate within the vector backbone; iii) the small T intron has been removed to prevent reduced reporter gene expression due to cryptic RNA splicing (22,23); and iv) a Kozak consensus sequence (25) has been inserted to increase the efficiency of translation initiation of the luciferase gene (17; Table 2).

Table 2. Changes Made to the pGL3 Vectors.

Change from pGL2	Purpose of Modification	Reference
Modifications made to the luciferase gene (<i>luc</i> to <i>luc+</i>).	Changes eliminate peroxisome targeting of expressed protein, eliminate consensus binding sequences for various genetic regulatory proteins, improve codon usage for mammalian and plant cells, and provide convenient restriction sites.	(18)
A unique <i>Nco</i> I site created at 5' end of <i>luc+</i> gene. <i>Nco</i> I sites removed from SV40 enhancer and promoter regions.	Ability to create N-terminal gene fusions with <i>luc+</i> using unique <i>Nco</i> I site.	
Intron from SV40 small T antigen removed.	Intron from SV40 small T antigen can reduce expression when placed 3' of certain genes due to cryptic splicing.	(23,24)
Poly(A) site for background reduction changed from SV40 early site to a synthetic poly(A) and transcriptional pause site.	Avoids possible recombination between two SV40 poly(A) sequences in the same plasmid.	(20,21)
Poly(A) signal for <i>luc+</i> changed from early to late SV40 poly(A) signal.	Late SV40 poly(A) signal is more efficient than early SV40 poly(A).	(19)
Kozak consensus sequence created immediately 5' of the <i>luc+</i> gene.	Provides for optimal translation efficiency.	(25)
Unique <i>Xba</i> I site created just downstream of the <i>luc+</i> gene.	User convenience; facilitates subcloning of the <i>luc+</i> gene.	
<i>Sma</i> I site moved to internal position in MCS.	User convenience; blunt-ended inserts can now be cleaved on either side by restriction endonucleases.	

D. Distinguishing Features of the pGL3 Luciferase Reporter Vectors

Maps of the pGL3-Basic, pGL3-Promoter, pGL3-Enhancer and pGL3-Control Vectors are shown in Figures 1–4. The DNA sequences and listings of restriction sites for these vectors are provided in Section X.H–K.

pGL3-Basic

The pGL3-Basic Vector lacks eukaryotic promoter and enhancer sequences, allowing maximum flexibility in cloning putative regulatory sequences. Expression of luciferase activity in cells transfected with this plasmid depends on insertion and proper orientation of a functional promoter upstream from *luc+*. Potential enhancer elements can also be inserted upstream of the promoter or in the *Bam*H I or *Sal* I sites downstream of the *luc+* gene.

pGL3-Enhancer

The pGL3-Enhancer Vector contains an SV40 enhancer located downstream of *luc+* and the poly(A) signal. This aids in the verification of functional promoter elements because the presence of an enhancer will often result in transcription of *luc+* at higher levels.

pGL3-Promoter

The pGL3-Promoter Vector contains an SV40 promoter upstream of the luciferase gene. DNA fragments containing putative enhancer elements can be inserted either upstream or downstream of the promoter-*luc+* transcriptional unit.

pGL3-Control

The pGL3-Control Vector contains SV40 promoter and enhancer sequences, resulting in strong expression of *luc+* in many types of mammalian cells. This plasmid is useful in monitoring transfection efficiency, in general, and is a convenient internal standard for promoter and enhancer activities expressed by pGL3 recombinants.

Note: The specific transcriptional characteristics of the pGL3 Vectors will vary for different cell types. This may be particularly true for COS cells, which contain the SV40 large T antigen. The SV40 large T antigen promotes replication from the SV40 origin, which is found in the promoter of the pGL3-Promoter and pGL3-Control Vectors. The combination of large T antigen and SV40 origin will result in a higher copy number of these vectors in COS cells, which in turn may result in increased expression of the reporter gene compared to other cell and vector combinations.

E. Mapping Genetic Elements Located Within DNA Fragments

The locations of functional elements within a DNA fragment are often determined by making a set of unidirectional nested deletions following the method of Henikoff (10) and then assaying for changes in biological activity. This method takes advantage of the unique properties of Exonuclease III (Exo III), which will digest 5' overhangs but not 3' overhangs or α -phosphorothioate nucleotide filled-in overhangs. Nested deletions of an insert DNA can be made directly in the pGL3 family of Reporter Vectors using this method, eliminating the need for subcloning steps. The multiple cloning site of the pGL3 Vectors contains

upstream *Kpn* I and *Sac* I restriction sites, which can be used to generate the 3' overhangs resistant to Exo III (Figures 1–5). After treatment with Exo III, S1 nuclease is added to remove the resulting ssDNA overhangs, and T4 DNA ligase is added to reclose the vectors. Deletion clones can be screened by gel electrophoresis of miniprep DNA, and the precise deletion endpoints within the promoter region can be determined by DNA sequencing using primers designed for the Luciferase Reporter Vectors.

F. Composition of Buffers and Solutions

lysis buffer for plasmid preps

25mM Tris-HCl (pH 7.5)
10mM EDTA
15% sucrose
2mg/ml lysozyme

M-9 plates (per liter)

6.0g Na₂HPO₄
3.0g KH₂PO₄
0.5g NaCl
1.0g NH₄Cl
15g agar

Add deionized water to approximately 1 liter. Autoclave. Cool to 50°C. Add the following sterilized solutions:

2.0ml 1M MgSO₄
0.1ml 1M CaCl₂
10.0ml 20% glucose (filter-sterilized)
1.0ml 1M thiamine-HCl

TE buffer (pH 8.0)

10mM Tris-HCl (pH 8.0)
1mM EDTA

TE-saturated

phenol:chloroform:isoamyl alcohol

Mix equal parts of TE buffer and phenol and allow the phases to separate. Then mix 1 part of the lower, phenol phase with 1 part of chloroform:isoamyl alcohol (24:1).

LB medium (per liter)

10g Bacto®-tryptone
5g Bacto®-yeast extract
5g NaCl

potassium acetate (pH 4.8)

Prepare 60ml of 5M potassium acetate. Add 11.5ml of glacial acetic acid and 28.5ml of water. This solution will be 3M with respect to potassium and 5M with respect to acetate. Store at 4°C.

G. References

1. *Protocols and Applications Guide*, Third Edition, (1996) Promega Corporation.
2. Sambrook, J. *et al.* (1989) *Molecular Cloning, A Laboratory Manual*, Cold Spring Harbor Press, Cold Spring Harbor, NY.
3. Brondyk, B. *et al.* (1994) A comparison of the Wizard™ Maxipreps DNA Purification System and alkaline lysis/cesium chloride method for isolating transfection-quality plasmid DNA. *Promega Notes* **47**, 2–5.
4. Birnboim, H.C. and Doly, J. (1979) A rapid alkaline extraction procedure for screening recombinant plasmid DNA. *Nucl. Acids Res.* **7**, 1513–23.
5. Schenborn, E. and Goiffon, V. (1991) Optimization of Transfectam®-mediated transfection using a luciferase reporter system. *Promega Notes* **33**, 8–11.
6. Cullen, B.R. (1987) Use of eukaryotic expression technology in the functional analysis of cloned genes. *Meth. Enzymol.* **152**, 684–704.
7. Ausubel, F.M. *et al.* (1988) *Current Protocols in Molecular Biology*, John Wiley and Sons, NY.
8. Rosenthal, N. (1987) Identification of regulatory elements of cloned genes with functional assays. *Meth. Enzymol.* **152**, 704–20.
9. Hawkins, E., Butler, B. and Wood, K.V. (2000) Bright-Glo™ and Steady-Glo™ Luciferase Assay Systems: Reagents for academic and industrial applications. *Promega Notes* **75**, 3–6.
10. Henikoff, S. (1987) Unidirectional digestion with exonuclease III in DNA sequence analysis. *Meth. Enzymol.* **155**, 156.
11. Dotto, G.P., Enea, V. and Zinder, N.D. (1981) Functional analysis of bacteriophage f1 intergenic region. *Virology* **114**, 463–73.
12. Dotto, G.P. and Zinder, N.D. (1983) The morphogenetic signal of bacteriophage f1. *Virology* **130**, 252–6.
13. Dotto, G.P., Huriuchi, K. and Zinder, N.D. (1984) The functional origin of bacteriophage f1 DNA replication. Its signals and domains. *J. Mol. Biol.* **172**, 507–21.
14. Hutchison, C.A. *et al.* (1978) Mutagenesis at a specific position in a DNA sequence. *J. Biol. Chem.* **253**, 6551–60.
15. Chien, A., Edgar, D.B. and Trela, J.M. (1976) Deoxyribonucleic acid polymerase from the extreme thermophile *Thermus aquaticus*. *J. Bacteriol.* **127**, 1550–7.
16. Kaledin, A.S., Sliusarenko, A.G. and Gorodetskii, S.T. (1980) Isolation and properties of DNA polymerase from extreme thermophilic bacteria *Thermus aquaticus* YT-1. *Biokhimiya* **45**, 644–51.
17. Innis, M.A. *et al.* (1988) DNA sequencing with *Thermus aquaticus* DNA polymerase and direct sequencing of polymerase chain reaction-amplified DNA. *Proc. Natl. Acad. Sci. USA* **85**, 9436–40.
18. Groskreutz, D.J. *et al.* (1995) Increased expression and convenience with the new pGL3 Luciferase Reporter Vectors. *Promega Notes* **50**, 2.
19. Sherf, B.A. and Wood, K.V. (1994) Firefly luciferase engineered for improved genetic reporting. *Promega Notes* **49**, 14–21.

20. Carswell, S. and Alwine, J.C. (1989) Efficiency of utilization of the simian virus 40 late polyadenylation site: Effects of upstream sequences. *Mol. Cell. Biol.* **9**, 4248–58.
21. Levitt, N. *et al.* (1989) Definition of an efficient synthetic poly(A) site. *Genes and Dev.* **3**, 1019–25.
22. Enriquez-Harris, P. *et al.* (1991) A pause site for RNA polymerase II is associated with termination of transcription. *EMBO J.* **10**, 1833–42.
23. Evans, M.J. and Scarpulla, R.C. (1989) Introns in the 3' untranslated region can inhibit chimeric CAT and beta-galactosidase gene expression. *Gene* **84**, 135–42.
24. Huang, M.T.F. and Gorman, C.M. (1990) The simian virus 40 small-t intron, present in many common expression vectors, leads to aberrant splicing. *Mol. Cell. Biol.* **10**, 1805–10.
25. Kozak, M. (1989) The scanning model for translation: An update. *J. Cell Biol.* **108**, 229–41.

H. pGL3-Basic Vector Restriction Sites and Sequence

The following restriction enzyme tables were constructed using DNASTAR® sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3' end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes, or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. Vector sequences are also available in the GenBank® database (GenBank®/EMBL Accession Number U47295) and on the Internet at www.promega.com/vectors/.

Table 3. Restriction Enzymes That Cut the pGL3-Basic Vector Between 1 and 5 Times.

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
Acc I	1	2011	Eag I	3	1755, 1759, 4651
Acc III	2	783, 1299	EclHK I	1	3153
Acc65 I	1	1	Eco47 III	1	2136
Acy I	4	95, 121, 1514, 3690	Eco52 I	3	1755, 1759, 4651
Afl III	3	15, 581, 2260	EcoCR I	1	9
Alw26 I	5	1111, 1343, 1409, 3214, 3990	EcoN I	3	645, 1045, 1705
Alw44 I	2	2574, 3820	Ehe I	1	122
AlwNI	1	2676	Fse I	1	1761
AspH I	5	11, 1553, 2578, 3739, 3824	Fsp I	2	3375, 4548
Ava I	3	26, 32, 1144	Hinc II	3	1392, 1902, 2012
Ava II	3	1267, 3291, 3513	Hind II	3	1392, 1902, 2012
BamH I	1	2004	Hind III	1	53
Ban II	4	11, 33, 1112, 4231	Hpa I	1	1902
Bbe I	1	124	Hsp92 I	4	95, 121, 1514, 3690
Bbs I	4	98, 1376, 1492, 2089	Kas I	1	120
Bbu I	1	751	Kpn I(p)	1	5
Bcl I	1	668	Mlu I	1	15
Bgl I	2	3273, 4541	Nae I	3	1759, 2130, 4199
Bgl II	1	36	Nar I	1	121
Bsa I	1	3214	Nco I	1	86
BsaA I	1	4302	NgoM IV	3	1757, 2128, 4197
BsaB I	1	2003	Nhe I	1	21
BsaH I	4	95, 121, 1514, 3690	Not I	1	4651
BsaM I	3	60, 1823, 1916	Nsp I	2	751, 2264
Bsm I	3	60, 1823, 1916	PaeR7 I	1	32
BspH I	3	671, 2980, 3988	PpuM I	1	1267
BspM I	3	1477, 1486, 4781	PshA I	1	2075
Bsr BR I	1	2003	Psp5 II	1	1267
BsrG I	1	578	PspA I	1	26
BssS I	2	2433, 3817	Pvu I	2	3523, 4569
BstZ I	3	1755, 1759, 4651	Sac I	1	11
Cla I	3	1997, 4709, 4813	Sal I	1	2010
Csp45 I	1	257	Sca I	3	253, 3633, 4716
Dra I	4	1963, 3019, 3038, 3730	SgrA I	1	1516
Dra II	1	1267	Sin I	3	1267, 3291, 3513
Dra III	1	4305	Sma I	1	28
Drd I	3	1489, 2368, 4349	Sph I	1	751
Dsa I	2	86, 458	Srf I	1	28
Eae I	4	1755, 1759, 3541, 4651	Ssp I	3	3957, 4510, 4625
			Sty I	1	86
			Vsp I	1	3325
			Xba I	1	1742
			Xcm I	1	823
			Xho I	1	32
			Xma I	1	26
			Xmn I	1	3752

Table 4. Restriction Enzymes That Do Not Cut the pGL3-Basic Vector.

<i>Aat</i> II	<i>Bbr</i> P I	<i>Bst</i>X I	<i>Nde</i> I	<i>Ppu</i> 10 I	<i>Spe</i> I
<i>Acc</i> B7 I	<i>Blp</i> I	<i>Bsu</i>36 I	<i>Nru</i> I	<i>Pst</i> I	<i>SpI</i>
<i>Afl</i> II	<i>Bpu</i> 1102 I	<i>Csp</i> I	<i>Nsi</i> I	<i>Pvu</i> II	<i>Sse</i> 8387 I
<i>Age</i> I	<i>Bsp</i> 120 I	<i>Eco</i> 72 I	<i>Pac</i> I	<i>Rsr</i> II	<i>Stu</i> I
<i>Apa</i> I	<i>Bss</i>H II	<i>Eco</i> 81 I	<i>Pfl</i> M I	<i>Sac</i> II	<i>Swa</i> I
<i>Asc</i> I	<i>Bst</i> 1107 I	<i>Eco</i>R I	<i>Pin</i> A I	<i>Sfi</i> I	<i>Tth</i>111 I
<i>Avr</i> II	<i>Bst</i>98 I	<i>Eco</i>R V	<i>Pme</i> I	<i>Sgf</i> I^(q)	
<i>Bal</i> I	<i>Bst</i>E II	<i>I-Ppo</i> I	<i>Pml</i> I	<i>Sna</i>B I	

Table 5. Restriction Enzymes That Cut the pGL3-Basic Vector 6 or More Times.

<i>Aci</i> I	<i>Bst</i>71 I	<i>Fnu</i> 4H I	<i>Hsp</i>92 II	<i>Msp</i>A1 I	<i>Scr</i> F I
<i>Alu</i> I	<i>Bst</i>O I	<i>Fok</i> I	<i>Mae</i> I	<i>Nci</i> I	<i>Sfa</i> N I
<i>Ban</i> I	<i>Bst</i> U I	<i>Hae</i> II	<i>Mae</i> II	<i>Nde</i> I	<i>Taq</i> I
<i>Bbv</i> I	<i>Cfo</i> I	<i>Hae</i> III	<i>Mae</i> III	<i>Nla</i> III	<i>Tfi</i> I
<i>Bsa</i>0 I	<i>Cfr</i> 10 I	<i>Hga</i> I	<i>Mbo</i> I	<i>Nla</i> IV	<i>Tru</i>9 I
<i>Bsa</i> J I	<i>Dde</i> I	<i>Hha</i> I	<i>Mbo</i> II	<i>Ple</i> I	<i>Xho</i> II
<i>Bsp</i>1286 I	<i>Dpn</i> I	<i>Hinf</i> I	<i>Mnl</i> I	<i>Rsa</i> I	
<i>Bsr</i> I	<i>Dpn</i> II	<i>Hpa</i> II	<i>Mse</i> I	<i>Sau</i>3A I	
<i>Bsr</i> S I	<i>Ear</i> I	<i>Hph</i> I	<i>Msp</i> I	<i>Sau</i>96 I	

Note: The enzymes listed in boldface type are available from Promega.

pGL3-Basic Vector Sequence

The strand shown is the same as the ssDNA strand produced by this vector and also corresponds to the mRNA synthesized from the *luc+* gene.

```

1  GGTACCGAGC TCTTACGCGT GCTAGCCCGG GCTCGAGATC TGCGATCTAA
51  GTAAGCTTGG CATTCCGGTA CTGTTGGTAA AGCCACCATG GAAGACGCCA
101 AAAACATAAA GAAAGGCCCG GCGCCATTCT ATCCGCTGGA AGATGGAACC
151 GCTGGAGAGC AACTGCATAA GGCTATGAAG AGATACGCCC TGGTTCCTGG
201 AACAAATTGCT TTTACAGATG CACATATCGA GGTGGACATC ACTTACGCTG
251 AGTACTTCGA AATGTCCGTT CGGTTGGCAG AAGCTATGAA ACGATATGGG
301 CTGAATACAA ATCACAGAAT CGTCGTATGC AGTGAAAAC TCTTTCAATT
351 CTTTATGCCG GTGTTGGGCG CGTTATTTAT CGGAGTTGCA GTTGCGCCCG
401 CGAACGACAT TTATAATGAA CGTGAATTGC TCAACAGTAT GGGCATTTCG
451 CAGCCTACCG TGGTGTTCGT TTCCAAAAAG GGGTTGCAA AAATTTTGAA
501 CGTGCAAAAA AAGCTCCCAA TCATCCAAAA AATTATTATC ATGGATTCTA
551 AAACGGATTA CCAGGGATTT CAGTCGATGT ACACGTTTCGT CACATCTCAT
601 CTACCTCCCG GTTTTAATGA ATACGATTTT GTGCCAGAGT CCTTCGATAG
651 GGACAAGACA ATTGCACTGA TCATGAACTC CTCTGGATCT ACTGGTCTGC
701 CTAAAGGTGT CGCTCTGCCT CATAGAACTG CCTGCGTGAG ATTCTCGCAT
751 GCCAGAGATC CTATTTTTGG CAATCAAATC ATTCCGGATA CTGCGATTTT

```

pGL3-Basic Vector Sequence (continued)

```

801  AAGTGTGTT  CCATTCCATC  ACGGTTTTGG  AATGTTTACT  ACACTCGGAT
851  ATTTGATATG  TGGATTTCTGA  GTCGTCTTAA  TGTATAGATT  TGAAGAAGAG
901  CTGTTTCTGA  GGAGCCTTCA  GGATTACAAG  ATTCAAAGTG  CGCTGCTGGT
951  GCCAACCCCTA  TTCTCCTTCT  TCGCCAAAAG  CACTCTGATT  GACAAATACG
1001 ATTTATCTAA  TTTACACGAA  ATTGCTTCTG  GTGGCGCTCC  CCTCTCTAAG
1051 GAAGTCGGGG  AAGCGGTTGC  CAAGAGGTTT  CATCTGCCAG  GTATCAGGCA
1101 AGGATATGGG  CTCACTGAGA  CTACATCAGC  TATTCTGATT  ACACCCGAGG
1151 GGGATGATAA  ACCGGGCGCG  GTCGGTAAAG  TTGTTCCATT  TTTTGAAGCG
1201 AAGTTGTGG  ATCTGGATAC  CGGAAAACG  CTGGGCGTTA  ATCAAAGAGG
1251 CGAACTGTGT  GTGAGAGGTC  CTATGATTAT  GTCCGGTTAT  GTAAACAATC
1301 CGGAAGCGAC  CAACGCCTTG  ATTGACAAGG  ATGGATGGCT  ACATTCTGGA
1351 GACATAGCTT  ACTGGGACGA  AGACGAACAC  TTCTTCATCG  TTGACCGCCT
1401 GAAGTCTCTG  ATTAAGTACA  AAGGCTATCA  GGTGGCTCCC  GCTGAATTGG
1451 AATCCATCTT  GCTCCAACAC  CCCAACATCT  TCGACGCAGG  TGTCGCAGGT
1501 CTTCCCGACG  ATGACGCCGG  TGAACTTCCC  GCCGCCGTTG  TTGTTTTGGA
1551 GCACGGAAAG  ACGATGACGG  AAAAAGAGAT  CGTGGATTAC  GTCGCCAGTC
1601 AAGTAACAAC  CGCGAAAAAG  TTGCGCGGAG  GAGTTGTGTT  TGTGGACGAA
1651 GTACCGAAAG  GTCTTACCGG  AAAACTCGAC  GCAAGAAAAA  TCAGAGAGAT
1701 CCTCATAAAG  GCCAAGAAGG  GCGGAAAGAT  CGCCGTGTAA  TTCTAGAGTC
1751 GGGGCGGCCG  GCCGCTTCGA  GCAGACATGA  TAAGATACAT  TGATGAGTTT
1801 GGACAAACCA  CAACTAGAAAT  GCAGTGAAAA  AAATGCTTTA  TTTGTGAAAT
1851 TTGTGATGCT  ATTGCTTTAT  TTGTAACCAT  TATAAGCTGC  AATAAACAAAG
1901 TTAACAACAA  CAATTGCATT  CATTTTATGT  TTCAGGTTCA  GGGGGAGGTG
1951 TGGGAGGTTT  TTTAAAGCAA  GTAAAACCTC  TACAAATGTG  GTAAAATCGA
2001 TAAGGATCCG  TCGACCGATG  CCCTTGAGAG  CCTTCAACCC  AGTCAGCTCC
2051 TTCCGGTGGG  CGCGGGGCAT  GACTATCGTC  GCCGCACTTA  TGA CTGTCTT
2101 CTTTATCATG  CAACTCGTAG  GACAGGTGCC  GGCAGCGCTC  TTCCGCTTCC
2151 TCGCTCACTG  ACTCGCTGCG  CTCGGTCGTT  CGGCTGCGGC  GAGCGGTATC
2201 AGCTCACTCA  AAGGCGGTAA  TACGGTTATC  CACAGAAATCA  GGGGATAACG
2251 CAGGAAAGAA  CATGTGAGCA  AAAGGCCAGC  AAAAGGCCAG  GAACCGTAAA
2301 AAGGCCGCGT  TGCTGGCGTT  TTTCCATAGG  CTCCGCCCCC  CTGACGAGCA
2351 TCACAAAAAT  CGACGCTCAA  GTCAGAGGTG  GCGAAACCCG  ACAGGACTAT
2401 AAAGATACCA  GGCCTTTCCC  CCTGGAAGCT  CCCTCGTGCG  CTCTCCTGTT
2451 CCGACCTGTC  CGCTTACCGG  ATACCTGTCC  GCCTTTCTCC  CTTCGGGAAG
2501 CGTGGCGCTT  TCTCATAGCT  CACGCTGTAG  GTATCTCAGT  TCGGTGTAGG
2551 TCGTTCGCTC  CAAGCTGGGC  TGTGTGCACG  AACCCCCCGT  TCAGCCCGAC

```

pGL3-Basic Vector Sequence (continued)

```

2601  CGCTGCGCCT  TATCCGGTAA  CTATCGTCTT  GAGTCCAACC  CGGTAAGACA
2651  CGACTTATCG  CCACTGGCAG  CAGCCACTGG  TAACAGGATT  AGCAGAGCGA
2701  GGTATGTAGG  CGGTGCTACA  GAGTTCTTGA  AGTGGTGGCC  TAACTACGGC
2751  TACACTAGAA  GAACAGTATT  TGGTATCTGC  GCTCTGCTGA  AGCCAGTTAC
2801  CTTTCGAAAA  AGAGTTGGTA  GCTCTTGATC  CGGCAAACAA  ACCACCGCTG
2851  GTAGCGGTGG  TTTTTTTGTT  TGCAAGCAGC  AGATTACGCG  CAGAAAAAAA
2901  GGATCTCAAG  AAGATCCTTT  GATCTTTTCT  ACGGGGTCTG  ACGCTCAGTG
2951  GAACGAAAAC  TCACGTTAAG  GGATTTTGGT  CATGAGATTA  TCAAAAAGGA
3001  TCTTCACCTA  GATCCTTTTA  AATTAAAAAT  GAAGTTTTAA  ATCAATCTAA
3051  AGTATATATG  AGTAAACTTG  GTCTGACAGT  TACCAATGCT  TAATCAGTGA
3101  GGCACCTATC  TCAGCGATCT  GTCTATTTTC  TTCATCCATA  GTTGCCTGAC
3151  TCCCCGTCGT  GTAGATAACT  ACGATACGGG  AGGGCTTACC  ATCTGGCCCC
3201  AGTGCTGCAA  TGATACCGCG  AGACCCACGC  TCACCGGCTC  CAGATTTATC
3251  AGCAATAAAC  CAGCCAGCCG  GAAGGGCCGA  GCGCAGAAGT  GGTCTGCAA
3301  CTTTATCCGC  CTCCATCCAG  TCTATTAATT  GTTGCCGGGA  AGCTAGAGTA
3351  AGTAGTTCGC  CAGTTAATAG  TTTGCGCAAC  GTTGTTGCCA  TTGCTACAGG
3401  CATCGTGGTG  TCACGCTCGT  CGTTTGGTAT  GGCTTCATTC  AGCTCCGGTT
3451  CCCAACGATC  AAGGCGAGTT  ACATGATCCC  CCATGTTGTG  CAAAAAAGCG
3501  GTTAGCTCCT  TCGGTCCTCC  GATCGTTGTC  AGAAGTAAGT  TGGCCGCAGT
3551  GTTATCACTC  ATGGTTATGG  CAGCACTGCA  TAATTCTCTT  ACTGTCATGC
3601  CATCCGTAAG  ATGCTTTTCT  GTGACTGGTG  AGTACTCAAC  CAAGTCATTC
3651  TGAGAATAGT  GTATGCGGCG  ACCGAGTTGC  TCTTGCCCGG  CGTCAATACG
3701  GGATAATACC  GCGCCACATA  GCAGAACTTT  AAAAGTGCTC  ATCATTTGAA
3751  AACGTTCTTC  GGGGCGAAAA  CTCTCAAGGA  TCTTACCGCT  GTTGAGATCC
3801  AGTTCGATGT  AACCACCTCG  TGCACCCAAC  TGATCTTCAG  CATCTTTTAC
3851  TTTCAACCAGC  GTTCTGGGT  GAGCAAAAAAC  AGGAAGGCAA  AATGCCGCAA
3901  AAAAGGGAAT  AAGGGCGACA  CGGAAATGTT  GAATACTCAT  ACTCTTCCTT
3951  TTTCAATATT  ATTGAAGCAT  TTATCAGGGT  TATTGTCTCA  TGAGCGGATA
4001  CATATTTGAA  TGATTTTAGA  AAAATAAACA  AATAGGGGTT  CCGGCACAT
4051  TTTCCCGAAA  AGTGCCACCT  GACGCGCCCT  GTAGCGGCGC  ATTAAGCGCG
4101  GCGGGTGTGG  TGGTTACGCG  CAGCGTGACC  GCTACACTTG  CCAGCGCCCT
4151  AGCGCCCGCT  CCTTTCGCTT  TCTTCCCTTC  CTTTCTCGCC  ACGTTCGCCG
4201  GCTTTCCCCG  TCAAGCTCTA  AATCGGGGGC  TCCCTTTAGG  GTTCCGATTT
4251  AGTGCTTTAC  GGCACCTCGA  CCCCAAAAAA  CTTGATTAGG  GTGATGGTTC
4301  ACGTAGTGGG  CCATCGCCCT  GATAGACGGT  TTTTCGCCCT  TTGACGTTGG
4351  AGTCCACGTT  CTTTAATAGT  GGACTCTTGT  TCCAAACTGG  AACAACTC
4401  AACCTATCT  CGGTCTATTC  TTTTGATTTA  TAAGGGATTT  TGCCGATTTT

```

pGL3-Basic Vector Sequence (continued)

```
4451  GGCCTATTGG  TTAAAAAATG  AGCTGATTTA  ACAAAAATTT  AACGCGAATT
4501  TTAACAAAAT  ATTAACGCTT  ACAATTTGCC  ATTCGCCATT  CAGGCTGCGC
4551  AACTGTTGGG  AAGGGCGATC  GGTGCGGGCC  TCTTCGCTAT  TACGCCAGCC
4601  CAAGCTACCA  TGATAAGTAA  GTAATATTAA  GGTACGGGAG  GTACTTGGAG
4651  CGGCCGCAAT  AAAATATCTT  TATTTTCATT  ACATCTGTGT  GTTGGTTTTT
4701  TGTGTGAATC  GATAGTACTA  ACATACGCTC  TCCATCAAAA  CAAAACGAAA
4751  CAAAACAAAC  TAGCAAATA  GGCTGTCCCC  AGTGCAAGTG  CAGGTGCCAG
4801  AACATTTCTC  TATCGATA
```

I. pGL3-Enhancer Vector Restriction Sites and Sequence

The following restriction enzyme tables were constructed using DNASTAR® sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3' end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes, or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. Vector sequences are also available in the GenBank® database (GenBank®/EMBL Accession Number U47297) and on the Internet at www.promega.com/vectors/.

Table 6. Restriction Enzymes that cut the pGL3-Enhancer Vector Between 1 and 5 Times.

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
Acc I	1	2257	Dsa I	2	86, 458
Acc III	2	783,1299	Eae I	4	1755, 1759, 3787, 4897
Acc65 I	1	1	Eag I	3	1755, 1759, 4897
Acy I	4	95, 121, 1514, 3936	Ec/HK I	1	3399
Afl III	3	15, 581, 2506	Eco47 III	1	2382
Alw26 I	5	1111, 1343, 1409, 3460, 4236	Eco52 I	3	1755, 1759, 4897
Alw44 I	2	2820, 4066	EcoCR I	1	9
AlwNI	1	2922	EcoNI	3	645, 1045, 1705
AspHI	5	11, 1553, 2824, 3985, 4070	Ehe I	1	122
Ava I	3	26, 32, 1144	Fse I	1	1761
Ava II	3	1267, 3537, 3759	Fsp I	2	3621, 4794
BamHI	1	2250	Hinc II	3	1392, 1902, 2258
Ban II	4	11, 33, 1112, 4477	Hind II	3	1392, 1902, 2258
Bbe I	1	124	Hind III	1	53
Bbs I	4	98, 1376, 1492, 2335	Hpa I	1	1902
Bbu I	3	751, 2108, 2180	Hsp92 I	4	95, 121, 1514, 3936
Bcl I	1	668	Kas I	1	120
Bgl I	2	3519, 4787	Kpn I	1	5
Bgl II	1	36	Mlu I	1	15
Bsa I	1	3460	Nae I	3	1759, 2376, 4445
BsaA I	1	4548	Nar I	1	121
BsaB I	1	2003	Nco I	1	86
BsaH I	4	95, 121, 1514, 3936	NgoM IV	3	1757, 2374, 4443
BsaM I	3	60, 1823, 1916	Nhe I	1	21
Bsm I	3	60, 1823, 1916	Not I	1	4897
BspHI	3	671, 3226, 4234	Nsi I	2	2106, 2178
BspMI	3	1477, 1486, 5027	Nsp I	4	751, 2108, 2180, 2510
Bsr BR I	1	2003	PaeR7 I	1	32
BsrGI	1	578	Ppu10 I	2	2102, 2174
BssSI	2	2679, 4063	PpuMI	1	1267
BstZ I	3	1755, 1759, 4897	PshA I	1	2321
Cla I	3	1997, 4955, 5059	Psp5 II	1	1267
Csp45 I	1	257	PspA I	1	26
Dra I	4	1963, 3265, 3284, 3976	Pvu I	2	3769, 4815
Dra II	1	1267	Sac I	1	11
Dra III	1	4551	Sal I	1	2256
Drd I	3	1489, 2614, 4595	Sca I	3	253, 3879, 4962
			SgrA I	1	1516
			Sin I	3	1267, 3537, 3759
			Sma I	1	28
			Sph I	3	751, 2108, 2180

Table 6. Restriction Enzymes That Cut the pGL3-Enhancer Vector Between 1 and 5 Times (continued).

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
<i>Srf</i> I	1	28	<i>Xba</i> I	1	1742
<i>Ssp</i> I	3	4203, 4756, 4871	<i>Xcm</i> I	1	823
<i>Sty</i> I	1	86	<i>Xho</i> I	1	32
<i>Vsp</i> I	1	3571	<i>Xma</i> I	1	26
			<i>Xmn</i> I	1	3998

Table 7. Restriction Enzymes That Do Not Cut the pGL3-Enhancer Vector.

<i>Aat</i> II	<i>Bbr</i> P I	<i>Bst</i> X I	<i>Nde</i> I	<i>Pvu</i> II	<i>Sse</i> 8387 I
<i>Acc</i> B7 I	<i>Blp</i> I	<i>Bsu</i>36 I	<i>Nru</i> I	<i>Rsr</i> II	<i>Stu</i> I
<i>Afl</i> II	<i>Bpu</i> 1102I	<i>Csp</i> I	<i>Pac</i> I	<i>Sac</i> II	<i>Swa</i> I
<i>Age</i> I	<i>Bsp</i> 120 I	<i>Eco</i> 72 I	<i>Pfl</i> M I	<i>Sfi</i> I	<i>Tth</i> III I
<i>Apa</i> I	<i>Bss</i>H II	<i>Eco</i> 81 I	<i>Pin</i> A I	<i>Sgf</i> I	
<i>Asc</i> I	<i>Bst</i> 1107 I	<i>Eco</i>R I	<i>Pme</i> I	<i>Sna</i>B I	
<i>Avr</i> II	<i>Bst</i>98 I	<i>Eco</i>R V	<i>Pml</i> I	<i>Spe</i> I	
<i>Bal</i> I	<i>Bst</i>E II	<i>I-Ppo</i> I	<i>Pst</i> I	<i>Spl</i> I	

Table 8. Restriction Enzymes That Cut the pGL3-Enhancer Vector 6 or More Times.

<i>Aci</i> I	<i>Bst</i>71 I	<i>Fnu</i> 4H I	<i>Hsp</i>92 II	<i>Msp</i>A1 I	<i>Scr</i> F I
<i>Alu</i> I	<i>Bst</i>O I	<i>Fok</i> I	<i>Mae</i> I	<i>Nci</i> I	<i>Sfa</i> N I
<i>Ban</i> I	<i>Bst</i> U I	<i>Hae</i> II	<i>Mae</i> II	<i>Nde</i> II	<i>Taq</i> I
<i>Bbv</i> I	<i>Cfo</i> I	<i>Hae</i> III	<i>Mae</i> III	<i>Nla</i> III	<i>Tfi</i> I
<i>Bsa</i>O I	<i>Cfr</i> 10 I	<i>Hga</i> I	<i>Mbo</i> I	<i>Nla</i> IV	<i>Tru</i>9 I
<i>Bsa</i> J I	<i>Dde</i> I	<i>Hha</i> I	<i>Mbo</i> II	<i>Ple</i> I	<i>Xho</i> II
<i>Bsp</i>1286 I	<i>Dpn</i> I	<i>Hinf</i> I	<i>Mnl</i> I	<i>Rsa</i> I	
<i>Bsr</i> I	<i>Dpn</i> II	<i>Hpa</i> II	<i>Mse</i> I	<i>Sau</i>3A I	
<i>Bsr</i>S I	<i>Ear</i> I	<i>Hph</i> I	<i>Msp</i> I	<i>Sau</i>96 I	

Note: The enzymes listed in boldface type are available from Promega.

pGL3-Enhancer Vector Sequence

The strand shown is the same as the ssDNA strand produced by this vector and also corresponds to the mRNA synthesized from the *luc+* gene.

```

1  GGTACCGAGC TCTTACGCGT GCTAGCCCGG GCTCGAGATC TGCGATCTAA
51  GTAAGCTTGG CATTCGCGTA CTGTTGGTAA AGCCACCATG GAAGACGCCA
101 AAAACATAAA GAAAGGCCCG GCGCCATTCT ATCCGCTGGA AGATGGAACC
151 GCTGGAGAGC AACTGCATAA GGCTATGAAG AGATACGCCC TGGTTCCTGG
201 AACAAATTGCT TTTACAGATG CACATATCGA GGTGGACATC ACTTACGCTG
251 AGTACTTCGA AATGTCCGTT CGGTTGGCAG AAGCTATGAA ACGATATGGG
301 CTGAATACAA ATCACAGAAT CGTCGTATGC AGTGAAAAC TCTTCAATT
351 CTTTATGCCG GTGTTGGGCG CGTTATTTAT CGGAGTTGCA GTTGCGCCCG
401 CGAACGACAT TTATAATGAA CGTGAATTGC TCAACAGTAT GGGCATTTCG
451 CAGCCTACCG TGGTGTTCGT TTCCAAAAG GGGTTGCAAA AAATTTTGAA

```



pGL3-Enhancer Vector Sequence (continued)

501 CGTGCAAAAA AAGCTCCCAA TCATCCAAAA AATTATTATC ATGGATTCTA
551 AAACGGATTA CCAGGGATTT CAGTCGATGT ACACGTTTCGT CACATCTCAT
601 CTACCTCCCG GTTTTAATGA ATACGATTTT GTGCCAGAGT CCTTCGATAG
651 GGACAAGACA ATTGCACTGA TCATGAACTC CTCTGGATCT ACTGGTCTGC
701 CTAAAGGTGT CGCTCTGCCT CATAGAACTG CCTGCGTGAG ATTCTCGCAT
751 GCCAGAGATC CTATTTTTGG CAATCAAATC ATTCCGGATA CTGCGATTTT
801 AAGTGTTGTT CCATTCCATC ACGGTTTTGG AATGTTTACT ACACTCGGAT
851 ATTTGATATG TGGATTTCGA GTCGTCTTAA TGTATAGATT TGAAGAAGAG
901 CTGTTTCTGA GGAGCCTTCA GGATTACAAG ATTCAAAGTG CGCTGCTGGT
951 GCCAACCTTA TTCTCCTTCT TCGCCTAAAAG CACTCTGATT GACAAATACG
1001 ATTTATCTAA TTTACACGAA ATTGCTTCTG GTGGCGCTCC CCTCTCTAAG
1051 GAAGTCGGGG AAGCGGTTGC CAAGAGGTTT CATCTGCCAG GTATCAGGCA
1101 AGGATATGGG CTCACTGAGA CTACATCAGC TATTCTGATT ACACCCGAGG
1151 GGGATGATAA ACCGGGCGCG GTCGGTAAAG TTGTTCCATT TTTTGAAGCG
1201 AAGGTTGTGG ATCTGGATAC CGGGAAAACG CTGGGCGTTA ATCAAAGAGG
1251 CGAACTGTGT GTGAGAGGTC CTATGATTAT GTCCGGTTAT GTAAACAATC
1301 CGGAAGCGAC CAACGCCCTG ATTGACAAGG ATGGATGGCT ACATTCTGGA
1351 GACATAGCTT ACTGGGACGA AGACGAACAC TTCTTCATCG TTGACCGCCT
1401 GAAGTCTCTG ATTAAGTACA AAGGCTATCA GGTGGCTCCC GCTGAATTGG
1451 AATCCATCTT GCTCCAACAC CCCAACATCT TCGACGCAGG TGTCGCAGGT
1501 CTTCCCGACG ATGACGCCGG TGAACCTCCC GCCGCCGTTG TTGTTTTGGA
1551 GCACGGAAAG ACGATGACGG AAAAAGAGAT CGTGGATTAC GTCGCCAGTC
1601 AAGTAACAAC CGCGAAAAG TTGCGCGGAG GAGTTGTGTT TGTGGACGAA
1651 GTACCGAAAG GTCTTACCGG AAAACTCGAC GCAAGAAAAA TCAGAGAGAT
1701 CCTCATAAAG GCCAAGAAGG GCGGAAAGAT CGCCGTGTAA TTCTAGAGTC
1751 GGGGCGGCCG GCCGCTTCGA GCAGACATGA TAAGATACAT TGATGAGTTT
1801 GGACAAACCA CAACTAGAAT GCAGTGAAAA AAATGCTTTA TTTGTGAAAT
1851 TTGTGATGCT ATTGCTTTAT TTGTAACCAT TATAAGCTGC AATAACAAG
1901 TTAACAACAA CAATTGCATT CATTTTATGT TTCAGGTTCA GGGGGAGGTG
1951 TGGGAGGTTT TTTAAAGCAA GTAAAACCTC TACAAATGTG GTAAAATCGA
2001 TAAGGATCTG AACGATGGAG CGGAGAATGG GCGGAACTGG GCGGAGTTAG
2051 GGGCGGGATG GGCGGAGTTA GGGGCGGGAC TATGGTTGCT GACTAATTGA
2101 GATGCATGCT TTGCATACTT CTGCCTGCTG GGGAGCCTGG GGACTTTCCA
2151 CACCTGGTTG CTGACTAATT GAGATGCATG CTTTGCATAC TTCTGCCTGC
2201 TGGGGAGCCT GGGGACTTTC CACACCCTAA CTGACACACA TTCCACAGCG
2251 GATCCGTCGA CCGATGCCCT TGAGAGCCTT CAACCCAGTC AGCTCCTTCC

pGL3-Enhancer Vector Sequence (continued)

```

2301  GGTGGGCGCG  GGGCATGACT  ATCGTCGCCG  CACTTATGAC  TGTCTTCTTT
2351  ATCATGCAAC  TCGTAGGACA  GGTGCCGCA  GCGCTCTTCC  GCTTCCTCGC
2401  TCACTGACTC  GCTGCGCTCG  GTCGTTCCGC  TCGGCGAGC  GGTATCAGCT
2451  CACTCAAAGG  CGGTAATACG  GTTATCCACA  GAATCAGGGG  ATAACGCAGG
2501  AAAGAACATG  TGAGCAAAAG  GCCAGCAAAA  GGCCAGGAAC  CGTAAAAAGG
2551  CCGCGTTGCT  GCGGTTTTTC  CATAGGCTCC  GCCCCCTGA  CGAGCATCAC
2601  AAAAATCGAC  GCTCAAGTCA  GAGGTGGCGA  AACCCGACAG  GACTATAAAG
2651  ATACCAGGCG  TTTCCCCTG  GAAGCTCCCT  CGTGCCTCT  CCTGTTCCGA
2701  CCCTGCCGCT  TACCGGATAC  CTGTCCGCT  TTCTCCCTTC  GGGAAGCGTG
2751  GCGCTTTCTC  ATAGCTCACG  CTGTAGGTAT  CTCAGTTCGG  TGTAGGTCGT
2801  TCGCTCCAAG  CTGGGCTGTG  TGCACGAACC  CCCCGTTCAG  CCCGACCGCT
2851  GCGCCTTATC  CGGTAACTAT  CGTCTTGAGT  CCAACCCGGT  AAGACACGAC
2901  TTATCGCCAC  TGGCAGCAGC  CACTGGTAAC  AGGATTAGCA  GAGCGAGGTA
2951  TGTAGGCGGT  GCTACAGAGT  TCTTGAAGTG  GTGGCCTAAC  TACGGCTACA
3001  CTAGAAGAAC  AGTATTTGGT  ATCTGCGCTC  TGCTGAAGCC  AGTTACCTTC
3051  GGAAAAAGAG  TTGGTAGCTC  TTGATCCGGC  AAACAAACCA  CCGCTGGTAG
3101  CGGTGGTTTT  TTTGTTTGCA  AGCAGCAGAT  TACGCGCAGA  AAAAAAGGAT
3151  CTCAAGAAGA  TCCTTTGATC  TTTTCTACGG  GGTCTGACGC  TCAGTGGAAC
3201  GAAACTCAC  GTTAAGGGAT  TTTGGTCATG  AGATTATCAA  AAAGGATCTT
3251  CACCTAGATC  CTTTTAAAT  AAAAATGAAG  TTTTAAATCA  ATCTAAAGTA
3301  TATATGAGTA  AACTTGGTCT  GACAGTTACC  AATGCTTAAT  CAGTGAGGCA
3351  CCTATCTCAG  CGATCTGTCT  ATTTCTGTTCA  TCCATAGTTG  CCTGACTCCC
3401  CGTCGTGTAG  ATAACTACGA  TACGGGAGGG  CTTACCATCT  GGCCCCAGTG
3451  CTGCAATGAT  ACCGCGAGAC  CCACGCTCAC  CGGCTCCAGA  TTTATCAGCA
3501  ATAAACCAGC  CAGCCGGAAG  GGCCGAGCGC  AGAAGTGGTC  CTGCAACTTT
3551  ATCCGCTCC  ATCCAGTCTA  TTAATTGTTG  CCGGAAGCT  AGAGTAAGTA
3601  GTTCGCCAGT  TAATAGTTTG  CGCAACGTTG  TTGCCATTGC  TACAGGCATC
3651  GTGGTGTCAC  GCTCGTCGTT  TGGTATGGCT  TCATTCAGCT  CCGTTCCCA
3701  ACGATCAAGG  CGAGTTACAT  GATCCCCCAT  GTTGTGCAA  AAAGCGGTTA
3751  GCTCCTTCGG  TCCTCCGATC  GTTGTGAGAA  GTAAGTTGGC  CGCAGTGTTA
3801  TCACTCATGG  TTATGGCAGC  ACTGCATAAT  TCTCTTACTG  TCATGCCATC
3851  CGTAAGATGC  TTTTCTGTGA  CTGGTGAGTA  CTCAACCAAG  TCATTCTGAG
3901  AATAGTGAT  GCGGCGACCG  AGTTGCTCTT  GCCCGGCGTC  AATACGGGAT
3951  AATACCGCGC  CACATAGCAG  AACTTTAAAA  GTGCTCATCA  TTGAAAACG
4001  TTCTTCGGGG  CGAAAACCTC  CAAGGATCTT  ACCGCTGTTG  AGATCCAGTT
4051  CGATGTAACC  CACTCGTGCA  CCCAACTGAT  CTTCAGCATC  TTTTACTTTC
4101  ACCAGCGTTT  CTGGGTGAGC  AAAAACAGGA  AGGCAAAATG  CCGCAAAAAA

```

pGL3-Enhancer Vector Sequence (continued)

```

4151  GGAATAAGG  GCGACACGGA  AATGTTGAAT  ACTCATACTC  TTCCTTTTTTC
4201  AATATTATTG  AAGCATTTAT  CAGGGTTATT  GTCTCATGAG  CGGATACATA
4251  TTTGAATGTA  TTTAGAAAAA  TAAACAAATA  GGGGTTCCGC  GCACATTTCC
4301  CCGAAAAGTG  CCACCTGACG  CGCCCTGTAG  CGGCGCATT  AGCGCGGCGG
4351  GTGTGGTGGT  TACGCGCAGC  GTGACCGCTA  CACTTGCCAG  CGCCCTAGCG
4401  CCCGCTCCTT  TCGCTTTCTT  CCCTTCCTTT  CTCGCCACGT  TCGCCGGCTT
4451  TCCCCGTCAA  GCTCTAAATC  GGGGGCTCCC  TTTAGGGTTC  CGATTTAGTG
4501  CTTTACGGCA  CCTCGACCCC  AAAAACTTG  ATTAGGGTGA  TGGTTCACGT
4551  AGTGGGCCAT  CGCCCTGATA  GACGGTTTTT  CGCCCTTTGA  CGTTGGAGTC
4601  CACGTTCTTT  AATAGTGGAC  TCTTGTTCCA  AACTGGAACA  AACTCAACC
4651  CTATCTCGGT  CTATTCTTTT  GATTTATAAG  GGATTTTGCC  GATTTCGGCC
4701  TATTGGTTAA  AAAATGAGCT  GATTTAACAA  AAATTTAACG  CGAATTTTAA
4751  CAAAATATTA  ACGCTTACAA  TTTGCCATT  GCCATTCAGG  CTGCGCAACT
4801  GTTGGGAAGG  GCGATCGGTG  CGGGCCTCTT  CGCTATTACG  CCAGCCCAAG
4851  CTACCATGAT  AAGTAAGTAA  TATTAAGGTA  CGGGAGGTAC  TTGGAGCGGC
4901  CGCAATAAAA  TATCTTTATT  TTCATTACAT  CTGTGTGTTG  GTTTTTTGTG
4951  TGAATCGATA  GACTAACAT  ACGCTCTCCA  TCAAAACAAA  ACGAAACAAA
5001  ACAAACTAGC  AAAATAGGCT  GTCCCCAGTG  CAAGTGCAGG  TGCCAGAACA
5051  TTTCTCTATC  GATA

```

J. pGL3-Promoter Vector Restriction Sites and Sequence

The following restriction enzyme tables were constructed using DNASTAR® sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3' end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes, or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. Vector sequences are also available in the GenBank® database (GenBank®/EMBL Accession Number U47298) and on the Internet at www.promega.com/vectors/.

Table 9. Restriction Enzymes That Cut the pGL3-Promoter Vector Between 1 and 5 Times.

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
Acc I	1	2203	Eae I	4	1947, 1951, 3733, 4843
Acc III	2	975, 1491	Eag I	3	1947, 1951, 4843
Acc65 I	1	1	EcHK I	1	3345
Acy I	4	28, 313, 1706, 3882	Eco47 III	1	2328
Afl III	3	15, 773, 2452	Eco52 I	3	1947, 1951, 4843
Alw26 I	5	1303, 1535, 1601, 3406, 4182	EcoICR I	1	9
Alw44 I	2	2766, 4012	EcoN I	3	837, 1237, 1897
AlwN I	1	2868	Ehe I	1	314
AspH I	5	11, 1745, 2770, 3931, 4016	Fse I	1	1953
Ava I	3	26, 32, 1336	Fsp I	2	3567, 4740
Ava II	3	1459, 3483, 3705	Hinc II	3	1584, 2094, 2204
Avr II	1	229	Hind II	3	1584, 2094, 2204
BamH I	1	2196	Hind III	1	245
Ban II	4	11, 33, 1304, 4423	Hpa I	1	2094
Bbe I	1	316	Hsp92 I	4	287, 313, 1706, 3882
Bbs I	4	290, 1568, 1684, 2281	Kas I	1	312
Bbu I	1	943	Kpn I	1	5
Bcl I	1	860	Mlu I	1	15
Bgl I	3	182, 3465, 4733	Nae I	3	1951, 2322, 4391
Bgl II	1	36	Nar I	1	313
Bsa I	1	3406	Nco I	1	278
BsaA I	1	4494	Ngom IV	3	1949, 2320, 4389
BsaB I	2	48, 2195	Nhe I	1	21
BsaH I	4	287, 313, 1706, 3882	Not I	1	4843
BsaM I	3	252, 2015, 2108	Nsp I	2	943, 2456
Bsm I	3	252, 2015, 2108	PaeR7 I	1	32
BspH I	3	863, 3172, 4180	PpuM I	1	1459
BspM I	3	1669, 1678, 4973	PshA I	1	2267
Bsr BR I	2	48, 2195	Psp5 II	1	1459
BsrG I	1	770	PspA I	1	26
BssS I	2	2625, 4009	Pvu I	2	3715, 4761
BstZ I	3	1947, 1951, 4843	Sac I	1	11
Cla I	3	2189, 4901, 5005	Sal I	1	2202
Csp45 I	1	449	Sca I	3	445, 3825, 4908
Dra I	4	2155, 3211, 3230, 3922	Sfi I	1	182
Dra II	1	1459	SgrA I	1	1708
Dra III	1	4497	Sin I	3	1459, 3483, 3705
Drd I	3	1681, 2560, 4541	Sma I	1	28
Dsa I	2	278, 650	Sph I	1	943
			Srf I	1	28
			Ssp I	3	4149, 4702, 4817
			Stu I	1	228

Table 9. Restriction Enzymes That Cut the pGL3-Promoter Vector Between 1 and 5 Times (continued).

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
Sty I	2	229, 278	Xho I	1	32
Vsp I	1	3517	Xma I	1	26
Xba I	1	1934	Xmn I	1	3944
Xcm I	1	1015			

Table 10. Restriction Enzymes That Do Not Cut the pGL3-Promoter Vector.

Aat II	<i>Blp I</i>	Bsu36 I	Nru I	Pst I	<i>Sse8387 I</i>
Acc B7 I	<i>Bpu 1102 I</i>	Csp I	Nsi I	Pvu II	<i>Swa I</i>
<i>Afl II</i>	<i>Bsp120 I</i>	<i>Eco72 I</i>	<i>Pac I</i>	<i>Rsr II</i>	Tth111 I
Age I	BssH II	<i>Eco81 I</i>	<i>PflM I</i>	Sac II	
Apa I	<i>Bst1 107 I</i>	EcoR I	<i>PinA I</i>	Sgf I	
<i>Asc I</i>	Bst 98 I	EcoR V	<i>Pme I</i>	SnaB I	
Bal I	Bst E II	I-Ppo I	<i>Pml I</i>	Spe I	
<i>BbrP I</i>	Bst X I	Nde I	<i>Ppu10 I</i>	<i>SpI</i>	

Table 11. Restriction Enzymes That Cut the pGL3-Promoter Vector 6 or More Times.

<i>Aci I</i>	Bst7 1 I	<i>Fnu4H I</i>	Hsp92 II	MspA1 I	<i>ScrF I</i>
Alu I	Bst O I	Fok I	<i>Mae I</i>	Nci I	<i>SfaN I</i>
Ban I	<i>BstU I</i>	Hae II	<i>Mae II</i>	Nde II	Taq I
<i>Bbv I</i>	Cfo I	Hae III	<i>Mae III</i>	<i>Nla III</i>	<i>Tfi I</i>
Bsa0 I	<i>Cfr10 I</i>	<i>Hga I</i>	Mbo I	<i>Nla IV</i>	Tru9 I
<i>BsaJ I</i>	Dde I	Hha I	Mbo II	<i>Ple I</i>	Xho II
Bsp1286 I	Dpn I	Hinf I	<i>Mnl I</i>	Rsa I	
<i>Bsr I</i>	<i>Dpn II</i>	Hpa II	<i>Mse I</i>	Sau3A I	
Bsr S I	<i>Ear I</i>	<i>Hph I</i>	Msp I	Sau96 I	

Note: The enzymes listed in boldface type are available from Promega.

pGL3-Promoter Vector Sequence

The strand shown is the same as the ssDNA strand produced by this vector and corresponds to the mRNA synthesized from the *luc+* gene.

```

1  GGTACCGAGC  TCTTACGCGT  GCTAGCCCGG  GCTCGAGATC  TGCGATCTGC
51  ATCTCAATTA  GTCAGCAACC  ATAGTCCCGC  CCCTAACTCC  GCCCATCCCG
101 CCCCTAACTC  CGCCCAGTTC  CGCCCATTCT  CCGCCCCATC  GCTGACTAAT
151 TTTTTTTTATT  TATGCAGAGG  CCGAGGCCGC  CTCGGCCTCT  GAGCTATTCC
201 AGAAGTAGTG  AGGAGGCTTT  TTTGGAGGCC  TAGGCTTTTG  CAAAAGCTT
251 GGCATTCCGG  TACTGTTGGT  AAAGCCACCA  TGAAGACGC  CAAAACATA
301 AAGAAAGGCC  CGGCGCCATT  CTATCCGCTG  GAAGATGGAA  CCGCTGGAGA
351 GCAACTGCAT  AAGGCTATGA  AGAGATACGC  CCTGGTTCCT  GGAACAATTG
401 CTTTTACAGA  TGCACATATC  GAGGTGGACA  TCACTTACGC  TGAGTACTTC
451 GAAATGTCCG  TTCGGTTGGC  AGAAGCTATG  AAACGATATG  GGCTGAATAC

```

pGL3-Promoter Vector Sequence (continued)

```

501  AAATCACAGA  ATCGTCGTAT  GCAGTGAAAA  CTCTCTTCAA  TTCTTTATGC
551  CGGTGTTGGG  CGCGTTATTT  ATCGGAGTTG  CAGTTGCGCC  CGCGAACGAC
601  ATTTATAATG  AACGTGAATT  GCTCAACAGT  ATGGGCATTT  CGCAGCCTAC
651  CGTGGTGTTT  GTTTCCAAAA  AGGGGTTGCA  AAAAAATTTG  AACGTGCAAA
701  AAAAGCTCCC  AATCATCCAA  AAAATTATTA  TCATGGATTG  TAAAACGGAT
751  TACCAGGGAT  TTCAGTCGAT  GTACACGTTT  GTCACATCTC  ATCTACCTCC
801  CGGTTTTAAT  GAATACGATT  TTGTGCCAGA  GTCCTTCGAT  AGGGACAAGA
851  CAATTGCACT  GATCATGAAC  TCCTCTGGAT  CTA CTGGTCT  GCCTAAAGGT
901  GTCGCTCTGC  CTCATAGAAC  TGCCTGCGTG  AGATTCTCGC  ATGCCAGAGA
951  TCCTATTTTT  GGCAATCAAA  TCATTCCGGA  TACTGCGATT  TTAAGTGTG
1001 TTCCATTCCA  TCACGGTTTT  GGAATGTTTA  CTACACTCGG  ATATTTGATA
1051 TGTGGATTTT  GAGTCGTCTT  AATGTATAGA  TTTGAAGAAG  AGCTGTTTCT
1101 GAGGAGCCTT  CAGGATTACA  AGATTCAAAG  TGCCTGCTGT  GTGCCAACCC
1151 TATTCTCCTT  CTTTCGCCAA  AGCACTCTGA  TTGACAAATA  CGATTTATCT
1201 AATTTACACG  AAATTGCTTC  TGGTGGCGCT  CCCCTCTCTA  AGGAAGTCGG
1251 GGAAGCGGTT  GCCAAGAGGT  TCCATCTGCC  AGGTATCAGG  CAAGGATATG
1301 GGCTCACTGA  GACTACATCA  GCTATTCTGA  TTACACCCGA  GGGGGATGAT
1351 AAACCGGGCG  CGGTCGGTAA  AGTTGTTCCA  TTTTTTGAAG  CGAAGGTTGT
1401 GGATCTGGAT  ACCGGGAAAA  CGCTGGGCGT  TAATCAAAGA  GGCGAACTGT
1451 GTGTGAGAGG  TCCTATGATT  ATGTCCGGTT  ATGTAACAA  TCCGGAAGCG
1501 ACCAACGCCT  TGATTGACAA  GGATGGATGG  CTACATTCTG  GAGACATAGC
1551 TTACTGGGAC  GAAGACGAAC  ACTTCTTCAT  CGTTGACCGC  CTGAAGTCTC
1601 TGATTAAGTA  CAAAGGCTAT  CAGGTGGCTC  CCGCTGAATT  GGAATCCATC
1651 TTGCTCCAAC  ACCCAACAT  CTTCGACGCA  GGTGTCGCAG  GTCTTCCCGA
1701 CGATGACGCC  GGTGAACTTC  CCGCCGCCGT  TGTTGTTTTG  GAGCACGGAA
1751 AGACGATGAC  GGAAAAGAG  ATCGTGGATT  ACGTCGCCAG  TCAAGTAACA
1801 ACCGCGAAAA  AGTTGCGCGG  AGGAGTTGTG  TTTGTGGACG  AAGTACCGAA
1851 AGGTCTTACC  GGAAAACCTG  ACGCAAGAAA  AATCAGAGAG  ATCCTCATAA
1901 AGGCCAAGAA  GGGCGGAAAG  ATCGCCGTGT  AATTCTAGAG  TCGGGGCGGC
1951 CGGCCGCTTC  GAGCAGACAT  GATAAGATAC  ATTGATGAGT  TTGGACAAAC
2001 CACAAC TAGA  ATGCAGTGAA  AAAAATGCTT  TATTTGTGAA  ATTTGTGATG
2051 CTATTGCTTT  ATTTGTAACC  ATTATAAGCT  GCAATAACA  AGTTAACAAAC
2101 AACAAATGCA  TTCATTTTAT  GTTTCAGGTT  CAGGGGGAGG  TGTGGGAGGT
2151 TTTTTAAAGC  AAGTAAAACC  TCTACAAATG  TGGTAAAATC  GATAAGGATC
2201 CGTCGACCGA  TGCCCTTGAG  AGCCTTCAAC  CCAGTCAGCT  CCTTCCGGTG
2251 GGCGCGGGGC  ATGACTATCG  TCGCCGCACT  TATGACTGTC  TTCTTTATCA
2301 TGCAACTCGT  AGGACAGGTG  CCGGCAGCGC  TCTTCCGCTT  CCTCGCTCAC

```



pGL3-Promoter Vector Sequence (continued)

2351 TGA^{CT}CGCTG CG^{CT}CGGTCG TTCGGCTGCG GCGAGCGGTA TCAGCTCACT
2401 CAAAGGCGGT AATACGGTTA TCCACAGAAAT CAGGGGATAA CGCAGGAAAG
2451 AACATGTGAG CAAAAGGCCA GCAAAGGCC AGGAACCGTA AAAAGGCCGC
2501 GTTGCTGGCG TTTTCCATA GGCTCCGCC CCCTGACGAG CATCACAAAA
2551 ATCGACGCTC AAGTCAGAGG TGGCGAAACC CGACAGGACT ATAAAGATAC
2601 CAGGCGTTTC CCCCTGGAAG CTCCCTCGTG CGCTCTCCTG TTCCGACCCT
2651 GCCGCTTACC GGATACCTGT CCGCCTTTCT CCCTTCGGGA AGCGTGGCGC
2701 TTTCTCATAG CTCACGCTGT AGGTATCTCA GTTCGGTGTA GGTCGTTTCG
2751 TCCAAGCTGG GCTGTGTGCA CGAACCCCC GTTCAGCCCG ACCGCTGCGC
2801 CTTATCCGGT AACTATCGTC TTGAGTCCAA CCCGGTAAGA CACGACTTAT
2851 CGCCACTGGC AGCAGCCACT GGTAACAGGA TTAGCAGAGC GAGGTATGTA
2901 GGCGGTGCTA CAGAGTTCTT GAAGTGGTGG CCTAACTACG GCTACACTAG
2951 AAGAACAGTA TTTGGTATCT GCGCTCTGCT GAAGCCAGTT ACCTTCGGAA
3001 AAAGAGTTGG TAGCTCTTGA TCCGGCAAAC AAACCACCGC TGGTAGCGGT
3051 GGTTTTTTTTG TTTGCAAGCA GCAGATTACG CGCAGAAAAA AAGGATCTCA
3101 AGAAGATCCT TTGATCTTTT CTACGGGGTC TGACGCTCAG TGGAACGAAA
3151 ACTCACGTTA AGGGATTTTG GTCATGAGAT TATCAAAAAG GATCTTCACC
3201 TAGATCCTTT TAAATTAATA ATGAAGTTTT AAATCAATCT AAAGTATATA
3251 TGAGTAAACT TGGTCTGACA GTTACCAATG CTTAATCAGT GAGGCACCTA
3301 TCTCAGCGAT CTGTCTATTT CGTTCATCCA TAGTTGCCTG ACTCCCCGTC
3351 GTGTAGATAA CTACGATACG GGAGGGCTTA CCATCTGGCC CCAGTGCTGC
3401 AATGATACCG CGAGACCCAC GCTCACCGGC TCCAGATTTA TCAGCAATAA
3451 ACCAGCCAGC CGGAAGGGCC GAGCGCAGAA GTGGTCCTGC AACTTTATCC
3501 GCCTCCATCC AGTCTATTAA TTGTTGCCGG GAAGCTAGAG TAAGTAGTTC
3551 GCCAGTTAAT AGTTTGCGCA ACGTTGTTGC CATTGCTACA GGCATCGTGG
3601 TGTCACGCTC GTCGTTTGGT ATGGCTTCAT TCAGCTCCGG TTTCCAACGA
3651 TCAAGGCGAG TTACATGATC CCCCATGTTG TGCAAAAAAG CGGTTAGCTC
3701 CTTTCGGTCCT CCGATCGTTG TCAGAAGTAA GTTGGCCGCA GTGTTATCAC
3751 TCATGGTTAT GGCAGCACTG CATAATTCTC TFACTGTCAT GCCATCCGTA
3801 AGATGCTTTT CTGTGACTGG TGAGTACTCA ACCAAGTCAT TCTGAGAATA
3851 GTGTATGCGG CGACCGAGTT GCTCTTGCCC GCGTCAATA CGGGATAATA
3901 CCGCGCCACA TAGCAGAACT TTAAAAGTGC TCATCATTTG AAAACGTTCT
3951 TCGGGGCGAA AACTCTCAAG GATCTTACCG CTGTTGAGAT CCAGTTCGAT
4001 GTAACCCACT CGTGACCCA ACTGATCTTC AGCATCTTTT ACTTTCACCA
4051 GCGTTTCTGG GTGAGCAAAA ACAGGAAGGC AAAATGCCGC AAAAAAGGGA
4101 ATAAGGGCGA CACGGAAATG TTGAATACTC AACTCTTCC TTTTCAATA
4151 TTATTGAAGC ATTTATCAGG GTTATTGTCT CATGAGCGGA TACATATTTG

pGL3-Promoter Vector Sequence (continued)

```

4201  AATGTATTTA  GAAAAATAAA  CAAATAGGGG  TTCCGCGCAC  ATTTCCCCGA
4251  AAAGTGCCAC  CTGACGCGCC  CTGTAGCGGC  GCATTAAGCG  CGGCGGGTGT
4301  GGTGGTTACG  CGCAGCGTGA  CCGCTACACT  TGCCAGCGCC  CTAGCGCCCCG
4351  CTCCTTTCGC  TTTCTTCCCT  TCCTTTCTCG  CCACGTTTCG  CGGCTTTCCC
4401  CGTCAAGCTC  TAAATCGGGG  GCTCCCTTTA  GGGTTCCGAT  TTAGTGCTTT
4451  ACGGCACCTC  GACCCCAAAA  AACTTGATTA  GGGTGATGGT  TCACGTAGTG
4501  GGCCATCGCC  CTGATAGACG  GTTTTTCGCC  CTTTGACGTT  GGAGTCCACG
4551  TTCTTTAATA  GTGGACTCTT  GTTCCAAACT  GGAACAACAC  TCAACCCTAT
4601  CTCGGTCTAT  TCTTTTGATT  TATAAGGGAT  TTTGCCGATT  TCGGCCTATT
4651  GGTTAAAAAA  TGAGCTGATT  TAACAAAAT  TTAACGCGAA  TTTTAACAAA
4701  ATATTAACGC  TTACAATTTG  CCATTCGCCA  TTCAGGCTGC  GCAACTGTTG
4751  GGAAGGGCGA  TCGGTGCGGG  CCTCTTCGCT  ATTACGCCAG  CCCAAGCTAC
4801  CATGATAAGT  AAGTAATATT  AAGGTACGGG  AGGTACTTGG  AGCGGCCGCA
4851  ATAAAATATC  TTTATTTTCA  TTACATCTGT  GTGTTGGTTT  TTTGTGTGAA
4901  TCGATAGTAC  TAACATACGC  TCTCCATCAA  AACAAAACGA  AACAAAACAA
4951  ACTAGCAAAA  TAGGCTGTCC  CCAGTGCAAG  TGCAGGTGCC  AGAACATTTT
5001  TCTATCGATA

```

K. pGL3-Control Vector Restriction Sites and Sequence

The following restriction enzyme tables were constructed using DNASTAR® sequence analysis software. Please note that this information has not been verified by restriction digestion with each enzyme listed. The location given specifies the 3' end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes, or if you identify a discrepancy, please contact your local Promega Branch or Distributor. Vector sequences are also available in the GenBank® database (GenBank®/EMBL Accession Number U47296) and on the Internet at www.promega.com/vectors/.

Table 12. Restriction Enzymes That Cut the pGL3-Control Vector Between 1 and 5 Times.

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
Acc I	1	2449	EcoICR I	1	9
Acc III	2	975, 1491	EcoN I	3	837, 1237, 1897
Acc65 I	1	1	Ehe I	1	314
Acy I	4	287, 313, 1706, 4128	Fse I	1	1953
Afl III	3	15, 773, 2698	Fsp I	2	3813, 4986
Alw26 I	5	1303, 1535, 1601, 3652, 4428	Hinc II	3	1584, 2094, 2450
Alw44 I	2	3012, 4258	Hind II	3	1584, 2094, 2450
AlwN I	1	3114	Hind III	1	245
AspH I	5	11, 1745, 3016, 4177, 4262	Hpa I	1	2094
Ava I	3	26, 32, 1336	Hsp92 I	4	287, 313, 1706, 4128
Ava II	3	1459, 3729, 3951	Kas I	1	312
Avr II	1	229	Kpn I	1	5
BamH I	1	2442	Mlu I	1	15
Ban II	4	11, 33, 1304, 4669	Nae I	3	1951, 2568, 4637
Bbe I	1	316	Nar I	1	313
Bbs I	4	290, 1568, 1684, 2527	Nco I	1	278
Bbu I	3	943, 2300, 2372	NgoM IV	3	1949, 2566, 4635
Bcl I	1	860	Nhe I	1	21
Bgl I	3	182, 3711, 4979	Not I	1	5089
Bgl II	1	36	Nsi I	2	2298, 2370
Bsa I	1	3652	Nsp I	4	943, 2300, 2372, 2702
BsaA I	1	4740	PaeR7 I	1	32
BsaB I	2	48, 2195	Ppu10 I	2	2294, 2366
BsaH I	4	287, 313, 1706, 4128	PpuM I	1	1459
BsaM I	3	252, 2015, 2108	PshA I	1	2513
Bsm I	3	252, 2015, 2108	Psp5 II	1	1459
BspH I	3	863, 3418, 4426	PspA I	1	26
BspM I	3	1669, 1678, 5219	Pvu I	2	3961, 5007
Bsr BR I	2	48, 2195	Sac I	1	11
BsrG I	1	770	Sal I	1	2448
BssS I	2	2871, 4255	Sca I	3	445, 4071, 5154
BstZ I	3	1947, 1951, 5089	Sfi I	1	182
Cla I	3	2189, 5147, 5251	SgrA I	1	1708
Csp45 I	1	449	Sin I	3	1459, 3729, 3951
Dra I	4	2155, 3457, 3476, 4168	Sma I	1	28
Dra II	1	1459	Sph I	3	943, 2300, 2372
Dra III	1	4743	Srf I	1	28
Drd I	3	1681, 2806, 4787	Ssp I	3	4395, 4948, 5063
Dsa I	2	278, 650	Stu I	1	228
Eag I	3	1947, 1951, 5089	Sty I	2	229, 278
EclHK I	1	3591	Vsp I	1	3763
Eco47 III	1	2574	Xba I	1	1934
Eae I	4	1947, 1951, 3979,	Xcm I	1	1015
Eco52 I	3	1947, 1951, 5089	Xho I	1	32
			Xma I	1	26
			Xmn I	1	4190

Table 13. Restriction Enzymes That Do Not Cut the pGL3-Control Vector.

Aat II	<i>Bbr</i> P I	Bst E II	Eco R V	<i>Pme</i> I	Sna B I
Acc B7 I	<i>Blp</i> I	Bst X I	I-Ppo I	<i>Pml</i> I	Spe I
<i>Afl</i> II	<i>Bpu</i> 1102I	Bsu 36 I	Nde I	Pst I	<i>Spl</i> I
Age I	<i>Bsp</i> 120 I	Csp I	Nru I	Pvu II	<i>Sse</i> 8387 I
Apa I	Bss H II	<i>Eco</i> 72 I	<i>Pac</i> I	<i>Rsr</i> II	<i>Swa</i> I
<i>Asc</i> I	<i>Bst</i> 1107 I	<i>Eco</i> 81 I	<i>Pfl</i> M I	Sac II	Tth 111 I
Bal I	Bst 98 I	Eco R I	<i>Pin</i> A I	Sgf I	

Table 14. Restriction Enzymes That Cut the pGL3-Control Vector 6 or More Times.

<i>Aci</i> I	Bst 71 I	Fok I	<i>Mae</i> I	Nci I	<i>Sfa</i> N I
Alu I	Bst O I	Hae II	<i>Mae</i> II	Nde II	Taq I
Ban I	<i>Bst</i> U I	Hae III	<i>Mae</i> III	<i>Nla</i> III	<i>Tfi</i> I
<i>Bbv</i> I	Cfo I	<i>Hga</i> I	Mbo I	<i>Nla</i> IV	Tru 9 I
Bsa 0 I	<i>Cfr</i> 10 I	Hha I	Mbo II	<i>Ple</i> I	Xho II
<i>Bsa</i> J I	Dde I	Hinf I	<i>Mnl</i> I	Rsa I	
Bsp 1286 I	Dpn I	Hpa II	<i>Mse</i> I	Sau 3A I	
<i>Bsr</i> I	<i>Ear</i> I	<i>Hph</i> I	Msp	Sau 96 I	
Bsr S I	<i>Fnu</i> 4H I	Hsp 92 II	Msp A1 I	<i>Scr</i> F I	

Note: The enzymes listed in boldface type are available from Promega.

pGL3-Control Vector Sequence

The strand shown is the same as the ssDNA strand produced by this vector and also corresponds to the mRNA synthesized from the *luc+* gene.

```

1  GGTACCGAGC TCTTACGCGT GCTAGCCCGG GCTCGAGATC TGCGATCTGC
51  ATCTCAATTA GTCAGCAACC ATAGTCCCGC CCCTAACTCC GCCCATCCCG
101 CCCCTAACTC CGCCCAGTTC CGCCCATTCT CCGCCCCATC GCTGACTAAT
151 TTTTTTTTATT TATGCAGAGG CCGAGGCCGC CTCGGCCTCT GAGCTATTCC
201 AGAAGTAGTG AGGAGGCTTT TTTGGAGGCC TAGGCTTTTG CAAAAAGCTT
251 GGCATTCCGG TACTGTTGGT AAAGCCACCA TGGAAAGACGC CAAAAACATA
301 AAGAAAGGCC CGGCGCCATT CTATCCGCTG GAAGATGGAA CCGCTGGAGA
351 GCAACTGCAT AAGGCTATGA AGAGATACGC CCTGGTTCCT GGAACAATTG
401 CTTTTACAGA TGCACATATC GAGGTGGACA TCACTTACGC TGAGTACTTC
451 GAAATGTCCG TTCGGTTGGC AGAAGCTATG AAACGATATG GGCTGAATAC
501 AAATCACAGA ATCGTCGTAT GCAGTGAAAA CTCTCTTCAA TTCTTTATGC
551 CGGTGTTGGG CGCGTTATTT ATCGGAGTTG CAGTTGCGCC CGCGAACGAC
601 ATTTATAATG AACGTGAATT GCTCAACAGT ATGGGCATTT CGCAGCCTAC

```

pGL3-Control Vector Sequence (continued)

```

651  CGTGGTGTTT  GTTCCAAAA  AGGGGTTGCA  AAAAATTTTG  AACGTGCAAA
701  AAAAGCTCCC  AATCATCCAA  AAAATTATTA  TCATGGATTTC  TAAAACGGAT
751  TACCAGGGAT  TTCAGTCGAT  GTACACGTTT  GTCACATCTC  ATCTACCTCC
801  CGGTTTTAAT  GAATACGATT  TTGTGCCAGA  GTCCTTCGAT  AGGGACAAGA
851  CAATTGCACT  GATCATGAAC  TCCTCTGGAT  CTA CTGGTCT  GCCTAAAGGT
901  GTCGCTCTGC  CTCATAGAAC  TGCTGCGTG  AGATTCTCGC  ATGCCAGAGA
951  TCCTATTTTT  GGCAATCAA  TCATTCCGGA  TACTGCGATT  TTAAGTGTTG
1001 TTCCATTCCA  TCACGGTTTT  GGAATGTTTA  CTACACTCGG  ATATTTGATA
1051 TGTGGATTTC  GAGTCGTCTT  AATGTATAGA  TTTGAAGAAG  AGCTGTTTCT
1101 GAGGAGCCTT  CAGGATTACA  AGATTCAAAG  TGCCTGCTG  GTGCCAACCC
1151 TATTCTCCTT  CTTCGCCAAA  AGCACTCTGA  TTGACAAATA  CGATTTATCT
1201 AATTTACACG  AAATTGCTTC  TGGTGGCGCT  CCCCTCTCTA  AGGAAGTCGG
1251 GGAAGCGGTT  GCCAAGAGGT  TCCATCTGCC  AGGTATCAGG  CAAGGATATG
1301 GGCTCACTGA  GACTACATCA  GCTATTCTGA  TTACACCCGA  GGGGGATGAT
1351 AAACCGGGCG  CGGTCCGTAA  AGTTGTTCCA  TTTTTTGAAG  CGAAGGTTGT
1401 GGATCTGGAT  ACCGGGAAAA  CGCTGGGCGT  TAATCAAAGA  GGCGAACTGT
1451 GTGTGAGAGG  TCCTATGATT  ATGTCCGGTT  ATGTAAACAA  TCCGGAAGCG
1501 ACCAACGCCT  TGATTGACAA  GGATGGATGG  CTACATTCTG  GAGACATAGC
1551 TTA CTGGGAC  GAAGACGAAC  ACTTCTTCAT  CGTTGACCGC  CTGAAGTCTC
1601 TGATTAAGTA  CAAAGGCTAT  CAGGTGGCTC  CCGCTGAATT  GGAATCCATC
1651 TTGCTCCAAC  ACCCCAACAT  CTTCGACGCA  GGTGTCGCAG  GTCTTCCCGA
1701 CGATGACGCC  GGTGAACTTC  CCGCCGCCGT  TGTGTTTTTG  GAGCACGGAA
1751 AGACGATGAC  GGAAAAAGAG  ATCGTGGATT  ACGTCGCCAG  TCAAGTAACA
1801 ACCGCGAAAA  AGTTGCGCGG  AGGAGTTGTG  TTTGTGGACG  AAGTACCGAA
1851 AGGTCTTACC  GGAAACTCG  ACGCAAGAAA  AATCAGAGAG  ATCCTCATAA
1901 AGGCCAAGAA  GGGCGGAAAG  ATCGCCGTGT  AATTCTAGAG  TCGGGGCGGC
1951 CGGCCGCTTC  GAGCAGACAT  GATAAGATAC  ATTGATGAGT  TTGGACAAAC
2001 CACA ACTAGA  ATGCAGTGAA  AAAAATGCTT  TATTTGTGAA  ATTTGTGATG
2051 CTATTGCTTT  ATTTGTAACC  ATTATAAGCT  GCAATAAACA  AGTTAACAAC
2101 AAC AATTGCA  TTCATTTTAT  GTTTCAGGT  CAGGGGGAGG  TGTGGGAGGT
2151 TTTTTAAAGC  AAGTAAAACC  TCTACAAATG  TGGTAAAATC  GATAAGGATC
2201 TGAACGATGG  AGCGGAGAAT  GGGCGGAACT  GGGCGGAGTT  AGGGGCGGGA
2251 TGGGCGGAGT  TAGGGGCGGG  ACTATGGTTG  CTGACTAATT  GAGATGCATG
2301 CTTTGCATAC  TTCTGCCTGC  TGGGGAGCCT  GGGGACTTTC  CACACCTGGT
2351 TGCTGACTAA  TTGAGATGCA  TGCTTTGCAT  ACTTCTGCCT  GCTGGGGAGC
2401 CTGGGGACTT  TCCACACCCT  AACTGACACA  CATTCCACAG  CGGATCCGTC
2451 GACCGATGCC  CTTGAGAGCC  TTCAACCCAG  TCAGCTCCTT  CCGGTGGGCG
2501 CGGGGCATGA  CTATCGTCGC  CGCACTTATG  ACTGTCTTCT  TTATCATGCA

```

pGL3-Control Vector Sequence (continued)

```

2551  ACTCGTAGGA  CAGGTGCCGG  CAGCGCTCTT  CCGCTTCCTC  GCTCACTGAC
2601  TCGCTGCGCT  CGGTCGTTTC  GCTGCGGCGA  GCGGTATCAG  CTCACTCAA
2651  GGCGGTAATA  CGGTTATCCA  CAGAATCAGG  GGATAACGCA  GGAAAGAACA
2701  TGTGAGCAAA  AGGCCAGCAA  AAGGCCAGGA  ACCGTAAAAA  GGCCGCGTTG
2751  CTGGCGTTTT  TCCATAGGCT  CCGCCCCCTT  GACGAGCATC  AAAAAAATCG
2801  ACGCTCAAGT  CAGAGGTGGC  GAAACCCGAC  AGGACTATAA  AGATAACCAGG
2851  CGTTTCCCCC  TGAAGCTCC  CTCGTGCGCT  CTCCTGTTCC  GACCCTGCCG
2901  CTTACCGGAT  ACCTGTCCGC  CTTTCTCCCT  TCGGGAAGCG  TGGCGCTTTC
2951  TCATAGCTCA  CGCTGTAGGT  ATCTCAGTTC  GGTGTAGGTC  GTTCGCTCCA
3001  AGCTGGGCTG  TGTGCACGAA  CCCCCGTTT  AGCCCGACCG  CTGCGCCTTA
3051  TCCGGTAACT  ATCGTCTTGA  GTCCAACCCG  GTAAGACACG  ACTTATCGCC
3101  ACTGGCAGCA  GCCACTGGTA  ACAGGATTAG  CAGAGCGAGG  TATGTAGGCG
3151  GTGCTACAGA  GTTCTTGAAG  TGGTGGCCTA  ACTACGGCTA  CACTAGAAGA
3201  ACAGTATTTG  GTATCTGCGC  TCTGCTGAAG  CCAGTTACCT  TCGGAAAAG
3251  AGTTGGTAGC  TCTTGATCCG  GCAAACAAAC  CACCGCTGGT  AGCGGTGGTT
3301  TTTTGTTTG  CAAGCAGCAG  ATTACGCGCA  GAAAAAAGG  ATCTCAAGAA
3351  GATCCTTTGA  TCTTTTCTAC  GGGGTCTGAC  GCTCAGTGGA  ACGAAAATC
3401  ACGTTAAGGG  ATTTTGGTCA  TGAGATTATC  AAAAAGGATC  TTCACCTAGA
3451  TCCTTTTAAA  TTAATAATGA  AGTTTTAAAT  CAATCTAAAG  TATATATGAG
3501  TAAACTTGGT  CTGACAGTTA  CCAATGCTTA  ATCAGTGAGG  CACCTATCTC
3551  AGCGATCTGT  CTATTTTCGT  CATCCATAGT  TGCCTGACTC  CCCGTCGTGT
3601  AGATAACTAC  GATACGGGAG  GGCTTACCAT  CTGGCCCCAG  TGCTGCAATG
3651  ATACCGCGAG  ACCCACGCTC  ACCGGCTCCA  GATTTATCAG  CAATAAACCA
3701  GCCAGCCGGA  AGGGCCGAGC  GCAGAAGTGG  TCCTGCAACT  TTATCCGCCT
3751  CCATCCAGTC  TATTAATTGT  TGCCGGGAAG  CTAGAGTAAG  TAGTTCGCCA
3801  GTTAATAGTT  TGCGCAACGT  TGTTGCCATT  GCTACAGGCA  TCGTGGTGTC
3851  ACGCTCGTCG  TTTGGTATGG  CTTCAATCAG  CTCCGGTTCC  CAACGATCAA
3901  GGCGAGTTAC  ATGATCCCC  ATGTTGTGCA  AAAAAGCGGT  TAGCTCCTTC
3951  GGTCCTCCGA  TCGTTGTCAG  AAGTAAGTTG  GCCGCAGTGT  TATCACTCAT
4001  GGTTATGGCA  GCACTGCATA  ATTCTCTTAC  TGTCATGCCA  TCCGTAAGAT
4051  GCTTTTCTGT  GACTGGTGAG  TACTCAACCA  AGTCATTCTG  AGAATAGTGT
4101  ATGCGGCGAC  CGAGTTGCTC  TTGCCCGGCG  TCAATACGGG  ATAATACCGC
4151  GCCACATAGC  AGAACTTTAA  AAGTGCTCAT  CATTGGAAAA  CGTTCTTCGG
4201  GGCGAAAAC  CTCAAGGATC  TTACCGCTGT  TGAGATCCAG  TTCGATGTAA
4251  CCCACTCGTG  CACCCAAC  ATCTTACGCA  TCTTTTACTT  TCACCAGCGT
4301  TTCTGGGTGA  GCAAAAACAG  GAAGGCAAAA  TGCCGCAAAA  AAGGGAATAA
4351  GGGCGACACG  GAAATGTTGA  AACTCATAAC  TCTTCCTTTT  TCAATATTAT
4401  TGAAGCATT  ATCAGGGTTA  TTGTCTCATG  AGCGGATACA  TATTTGAATG

```

pGL3-Control Vector Sequence (continued)

```
4451 TATTTAGAAA AATAAACAAA TAGGGGTTCC GCGCACATTT CCCCGAAAAG
4501 TGCCACCTGA CGCGCCCTGT AGCGGCGCAT TAAGCGCGGC GGGTGTGGTG
4551 GTTACGCGCA GCGTGACCGC TACACTTGCC AGCGCCCTAG CGCCCGCTCC
4601 TTTTCGCTTTC TTCCCTTCCT TTCTCGCCAC GTTCGCCGGC TTTCCCCGTC
4651 AAGCTCTAAA TCGGGGGCTC CCTTTAGGGT TCCGATTTAG TGCTTTACGG
4701 CACCTCGACC CCAAAAAACT TGATTAGGGT GATGGTTCAC GTAGTGGGCC
4751 ATCGCCCTGA TAGACGGTTT TTCGCCCTTT GACGTTGGAG TCCACGTTCT
4801 TTAATAGTGG ACTCTTGTTT CAAACTGGAA CAACACTCAA CCCTATCTCG
4851 GTCTATTCTT TTGATTTATA AGGGATTTTG CCGATTTTCGG CCTATTGGTT
4901 AAAAAATGAG CTGATTTAAC AAAAATTTAA CGCGAATTTT AACAAAATAT
4951 TAACGCTTAC AATTTGCCAT TCGCCATTCA GGCTGCGCAA CTGTTGGGAA
5001 GGGCGATCGG TCGGGGCCTC TTCGCTATTA CGCCAGCCCA AGCTACCATG
5051 ATAAGTAAGT AATATTAAGG TACGGGAGGT ACTTGAGCG GCCGCAATAA
5101 AATATCTTTA TTTTCATTAC ATCTGTGTGT TGGTTTTTTG TGTGAATCGA
5151 TAGTACTAAC ATACGCTCTC CATCAAACA AAACGAAACA AAACAAACTA
5201 GCAAAATAGG CTGTCCCCAG TGCAAGTGCA GGTGCCAGAA CATTCTCTA
5251 TCGATA
```

- (a) U.S. Pat. No. 5,670,356 has been issued to Promega Corporation for a modified luciferase technology.
- (b) The method of recombinant expression of *Coleoptera* luciferase is covered by U.S. Pat. Nos. 5,583,024, 5,674,713 and 5,700,673.
- (c) The PCR process is covered by patents issued and applicable in certain countries. Promega does not encourage or support the unauthorized or unlicensed use of the PCR process.
- (d) Licensed under U.S. Pat. No. 5,075,430.
- (e) U.S. Pat. No. 5,981,235 and Australian Pat. No. 729932 have been issued to Promega Corporation for methods for isolating nucleic acids using alkaline protease. Other patents are pending. Australian Pat. No. 730718 has been issued to Promega Corporation for an improved filtration system and method. Other patents are pending.
- (f) U.S. Pat. Nos. 5,658,548 and 5,808,041 and Australian Pat. No. 689815 have been issued to Promega Corporation for nucleic acid purification on silica gel and glass mixtures. Other patents are pending.
- (g) U.S. Pat. No. 6,027,945 and Australian Pat. No. 732756 have been issued to Promega Corporation for methods of isolating biological target materials using silica magnetic particles. Other patents are pending. U.S. Pat. No. 6,194,562 has been issued to Promega Corporation for endotoxin reduction in nucleic acid purification. Other patents are pending.
- (h) Transfectam is a registered trademark of Promega Corporation, which sells the Transfectam[®] product for research purposes only under agreement with CNRS-Paris, owner by assignment of U. S. Pat. No. 5,171,678. The Transfectam[®] product was developed by J.P. Behr and J.P. Loeffler and is covered by the aforementioned patent.
- (i) The cationic lipid component of the TransFast[™] Transfection Reagent is covered by U.S. Pat. Nos. 5,824,812, 5,869,715 and pending foreign patents.
- (j) The cationic lipid component of the Tfx[™] Reagents is covered by U.S. Pat. Nos. 5,527,928, 5,744,625 and 5,892,071, Australian Pat. No. 704189 and other pending foreign patents.
- (k) U.S. Pat. Nos. 5,283,179, 5,641,641, 5,650,289, 5,814,471, Australian Pat. No. 649289 and European Pat. No. 0 553 234 have been issued to Promega Corporation for a firefly luciferase assay method, which affords greater light output with improved kinetics as compared to the conventional assay. Other patents are pending. Certain applications of this product may require licenses from others.
- (l) U.S. Pat. No. 5,955,363 has been issued to Promega Corporation for a vector for in vitro mutagenesis and use thereof.
- (m) The 7-deaza-dGTP component is licensed from Boehringer Mannheim GmbH under U.S. Pat. Nos. 4,804,748 and 5,480,980.
- (n) U.S. Pat. Nos. 5,523,206 and 5,654,149, Australian Pat. No. 671820 and European Pat. No. 0 610 615 B1 have been issued to Promega Corporation for non-radioactive DNA sequencing.
- (o) U.S. Pat. No. 5,108,892 has been issued to Promega Corporation for the use of a modified *Taq* DNA polymerase to determine DNA sequence and amplify DNA sequence.
- (p) Notice to Purchaser: Limited Use License
This product is sold under licensing arrangements between Promega Corporation and Invitrogen Corporation. The purchase price of this product includes limited, nontransferable rights under U.S. Pat. Nos. 5,082,784 and 5,192,675 owned by Invitrogen Corporation to use the product only for the internal research purposes of the purchaser. For information on purchasing a license to use the purchased product for purposes other than the internal research of the purchaser, contact the Director of Licensing, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008.
- (q) U.S. Pat. No. 5,391,487 has been issued to Promega Corporation for Restriction Endonuclease *Sgf I*.

© 1994–2002 Promega Corporation. All Rights Reserved.

Altered Sites, Steady-Glo, Erase-a-Base, *fmoI*, ProFection, Transfectam and Wizard are trademarks of Promega Corporation and are registered with the U.S. Patent and Trademark Office. Bright-Glo, SILVER SEQUENCE, Tfx and TransFast are trademarks of Promega Corporation.

Bacto is a registered trademark of Difco Laboratories, Detroit, Michigan. DNASTAR is a registered trademark of DNASTAR, Inc. GenBank is a registered trademark of the U.S. Dept. of Health and Human Services.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.

pGL3 Luciferase Reporter Vectors: *Experienced User's Protocol*

This quick protocol is intended as an easy-to-follow reminder for experienced users. Please follow the complete protocol (Sections IV.A–D) the first time you use the pGL3 Luciferase Reporter Vectors.

<p>Isolation of Plasmid DNA (Section IV.D)</p>	<ol style="list-style-type: none"> 1. Prepare an overnight culture in 250ml of LB medium containing 100µg/ml ampicillin. 2. Collect cells by centrifuging at 5,000 × <i>g</i> for 15 minutes at 4°C. Discard the supernatant. 3. Resuspend cells in 6ml of freshly prepared ice-cold lysis buffer. Mix carefully using a 10ml pipet. Incubate the tube containing the cells and lysis buffer in ice water for 20 minutes. 4. Add 12ml of 0.1N NaOH, 1% SDS (prepared fresh). Mix carefully and thoroughly by inversion. Do not vortex. 5. Add 7.5ml of potassium acetate solution (pH 4.8). Mix carefully by inversion. Incubate tube in ice water for 10 minutes. 6. Centrifuge at 12,000 × <i>g</i> for 15 minutes. Transfer the supernatant to a new tube, avoiding the white precipitate. Add 50µl of RNase A (1mg/ml stock) to the supernatant. Incubate for 20 minutes at 37°C. 7. Extract with one volume of TE-saturated phenol:chloroform:isoamyl alcohol. Centrifuge at 12,000 × <i>g</i> for 10 minutes. 8. Save the upper, aqueous phase and repeat the TE-saturated phenol:chloroform:isoamyl alcohol extraction described in Step 7. 9. Extract with one volume of chloroform:isoamyl alcohol (24:1) by vortexing for 1 minute. Centrifuge at 12,000 × <i>g</i> for 10 minutes. 10. Transfer the upper aqueous phase to a new tube and add 2 volumes of 100% ethanol. Centrifuge at 12,000 × <i>g</i> for 20 minutes. 11. Optional: Dissolve the pellet in 1.6ml of water. Add 0.4ml of 4M NaCl and mix. Add 2ml of 13% (w/v) polyethylene glycol and mix. Incubate the tube in ice water for 60 minutes. Centrifuge at 12,000 × <i>g</i> for 10 minutes. 12. Remove the supernatant and wash the pellet with 70% ethanol. Centrifuge at 12,000 × <i>g</i> for 5 minutes. 13. Dry the pellet under vacuum. Dissolve the pellet in water or TE buffer (100–500µl).
---	---